CLASSIC SHEATH™

Splittable Hemostatic Introducer System

en	Instructions for Use	2
cz	Návod k použití	4
da	Brugervejledning	6
de	Gebrauchsanweisung	8
el	Οδηγίες Χρήσης	10
es	Instrucciones de uso	12
et	Kasutusjuhend	14
fi	Käyttöohjeet	16
fr	Instructions d'utilisation	18
hu	Használati utasítás	20
it	Istruzioni per l'uso	22
jp	取扱説明書	24
ko	사용 지침	26
lt	Naudojimo instrukcijos	28
nl	Gebruiksaanwijzing	30
no	Bruksinstruksjoner	32
po	Instrukcja użycia	34
pt	Instruções de Utilização	36
sc	使用说明	38
sk	Pokyny na používanie	40
sv	Bruksanvisning	42
tu	Kullanım Talimatları	44

CLASSIC SHEATH™

Splittable Hemostatic Introducer System

Instructions for Use

**This device is intended for one time use only.
Read instructions prior to use.**

**CAUTION: Federal (U.S.A.) Law restricts this
device to sale by or on the order of a physician.**

For U.S.-California Only.

Proposition 65, a State of California voter initiative, requires the following notice:

WARNING: This product and its packaging have been sterilized with ethylene oxide. This packaging may expose you to ethylene oxide, a chemical known to the state of California to cause cancer or birth defects or other reproductive harm.

Indications

For the introduction of various types of pacing leads and catheters.

Warnings

- This product is sensitive to light. Do not use if stored outside the protective outer carton.
- Store in a cool, dark, and dry place.
- Infusion through the sideport can be done only after all air is removed from the unit.

Precautions

- Do not alter this device in any way.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- Aspiration and saline flushing of the sheath and dilator should be performed to help minimize the potential for air embolism and clot formation.
- Indwelling introducer sheaths should be internally supported by a catheter, pacing lead, or dilator.
- Dilators, catheters, and pacing leads should be removed slowly from the sheath. Rapid removal may damage the valve resulting in flow of blood or air through the valve.
- Never advance or withdraw the guidewire or sheath when resistance is met. Determine the cause by fluoroscopy and take remedial action.
- When injecting or aspirating through the sheath, use the sideport only.

Adverse Events

Adverse events may include, but are not limited to the following:

- | | |
|-----------------|----------------|
| • Air embolus | • Pneumothorax |
| • Blood loss | • Hemothorax |
| • Vessel damage | • Pacing lead |
| • Infection | displacement |

Use Sterile Technique A suggested procedure:

1. Peel open package and place contents in sterile field.
2. Prep skin and drape in area of anticipated veni-puncture.
3. Flush the needle, sheath, and dilator prior to use.
4. Insert the dilator into the sheath until the dilator cap attaches to the sheath hub.
5. Insert the needle into vessel. Verify the needle position by observing venous blood return.
6. Aspirate the needle using the syringe.
7. Remove the syringe and insert soft tip of the guidewire through the needle into the vessel. Advance the guidewire to required depth. Leave an appropriate amount of the guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire's entrance into the superior vena cava and right atrium is suggested.
8. Hold the guidewire in place and remove needle. Do not withdraw the guidewire back into the needle as this may result in separation of the guidewire.
9. Thread the dilator/sheath assembly over the guidewire.
10. Advance the dilator/sheath assembly with a twisting motion over the guidewire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.
11. Once assembly is fully introduced into the venous system, remove the dilator cap from the sheath by rocking the dilator cap off the sheath hub.
12. Slowly retract the guidewire and dilator, leaving the sheath in position. The sheath will reduce the loss of blood and the inadvertent aspiration of air.
13. Aspirate all air from the sheath by using a syringe connected to the sideport.
14. Flush the introducer with saline through the sideport. If the introducer is to remain in place during lead positioning and testing, periodic flushing is advised.
15. Introduce pacemaker lead or catheter through the sheath and advance it into position.
16. Flush sheath with saline immediately before peeling sheath away in order to minimize backbleeding.
17. Split sheath by sharply snapping the sheath hub tabs and peeling sheath tube apart while withdrawing it from the vessel.

TVI Instructions for Use

The Transvalvular Insertion Tool (TVI) is an optional device for use with the

CLASSICsheath systems

This device is intended for one time use only.
Read instructions prior to use.

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Indications

For the introduction of various types of pacing or defibrillator leads and catheters.

Warnings

Improper use of the TVI can cause air embolism and back bleeding.

Precautions

- Do not alter this device in any way.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- When using the TVI tool, lead size may not exceed the size of the TVI tool.
- When using the TVI, always keep the exposed proximal end covered to prevent air embolism and back bleeding.
- When retaining a guidewire, the TVI tool must be placed over the guidewire and not along the side of it.

Use Sterile Technique A suggested procedure:

1. A TVI is to be used at the physician's discretion to open the sheath valve for ease of lead placement.
2. CAUTION! When the TVI is inserted into the ClassicSheath valve housing, all hemostasis is lost and the risk of air embolization and back bleeding exists. Always keep the proximal exposed end of the TVI covered with your thumb while the TVI is in use.
3. To use the TVI, insert the distal end of the TVI into the valve housing by gently pushing the TVI into the sheath.
4. Hold your thumb over the proximal exposed opening of the TVI to prevent air embolism or back bleeding.
5. Advance the pacing lead or catheter through the TVI and into the sheath.
6. As soon as the pacing lead or catheter is resting inside the sheath, pull the TVI back out of the sheath valve housing.
7. The TVI may then either be peeled away or it may be temporarily left resting on the shaft of the pacing lead or catheter.
8. Once the TVI has been withdrawn from the sheath, aspirate the sheath through the sideport, until any air which may have entered the sheath during the procedure is removed. Flush with Heparinized saline.

Graphic Symbols for Medical Device Labeling

Single Use Only /
Do Not Reuse

Do not use if package is
opened or damaged

Sterilized Using
Ethylene Oxide

Store in a cool, dark and
dry place

Caution, read instructions
prior to use

Size

Contents (Numeral repre-sents quantity of units inside)

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Merit Medical Systems, Inc.
U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd.