# HeartSpan®

## Transseptal Needle

---

Transseptal Needle and Stylet Set  
Transseptální jehla se zavaděčem  
Transseptalt kanyleudstyr  
Transseptales Nadel- und Mandrin-Set  
Σετ διαδιαφραγματικής βελόνας και στυλεού  
Conjunto de aguja y estilete  
Transseptaalinen neula- ja mandriinisarja  
Ensemble stylet et aiguille transseptale  
Transzszeptális tű és szonda készlet  
Ago transettale e set stiletto  
経中隔針および探り針セット  
Transseptale naald en stiletset  
Transseptal kanyle- og sondesett  
Zestaw igły przezprzegrodowej i mandrynu  
Kit com Agulha Transeptal e Estilete  
Transseptal nål och sondset  
Transseptal İğne ve Stile Seti  
Трансептальная игла со стилетом  
Transseptálna ihla so styletom

en	Directions for Use	2
cs	Návod k použití	4
da	Brugsanvisning	6
de	Gebrauchsanweisung	8
el	Οδηγίες χρήσης	10
es	Instrucciones de uso	12
fi	Käyttöohjeet	14
fr	Instructions d'utilisation	16
hu	Használati útmutató	18
it	Istruzioni per l'uso	20
ja	使用説明書	22
nl	Richtlijnen voor gebruik	24
no	Bruksanvisning	26
pl	Sposób użycia	28
pt	Instruções de Utilização	30
sv	Bruksanvisning	32
tr	Kullanım Yönergeleri	34
ru	Указания по применению	36
sk	Návod na použitie	38

# en - Directions for Use

Read instructions for use in its entirety prior to use.

Single use only disposable medical device.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

## Device Description

The transseptal needle kit consists of an outer needle cannula and an inner stylet. The needle is comprised of flexible thin walled tubing with an ergonomic hub and stopcock attached to the proximal end. The stylet consists of a solid wire that when inserted in the needle protrudes beyond the distal tip of the cannula.

## Indications for Use

HeartSpan transseptal needle is used to create the primary puncture in the interatrial septum when passing an introducer and/or catheter through the septum from the right side of the heart to the left side.

## Contraindications

- Left atrial thrombus or tumor
- Inability to lie flat
- Substantial deformity of the spine or chest
- Distorted anatomy due to congenital heart disease
- Previous intra-septum patch

## Warnings

- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death. This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- Only those physicians who specialize in the practice of invasive cardiology techniques should use this device. The device should be restricted for use by those specialists trained to perform transseptal procedures.
- Maintain continuous pressure monitoring and repeated biplane fluoroscopy during positioning.
- Caution should be used in patients with small left atrium, to avoid left atrial wall puncture.
- The transseptal needle should never be advanced until the catheter is positioned correctly on the atrial septum.
- Always ensure that the transseptal needle has clearly entered the left atrial cavity by confirming distinct left atrial pressure and fluoroscopy of the needle tip before advancing the dilator, sheath or catheter.
- Do not remove a dilator, sheath or catheter that has been inadvertently advanced into the pericardial space until the patient is in surgery.
- Do not reuse this device. The device must be discarded after one use using acceptable medical practices and applicable local, state and federal laws and regulations.
- Follow instructions accompanying the transseptal introducers.

## Precautions

- Store in a cool, dark and dry place.
- Inspect all components prior to use.
- If resistance is met while advancing or withdrawing the introducer or guidewire, fluoroscopy should be used to determine the cause.
- Prior to use, ensure the appropriate catheter is being used with the transseptal needle.
- Administer appropriate levels of peri-procedural anticoagulation therapy for patients undergoing left-sided and transseptal cardiac procedures. Administer anticoagulation therapy during and post-procedure according to the institution's standards to minimize bleeding and thrombotic complications
- Prior to the procedure, the patient must be hemodynamically stable. Certain conditions may require special consideration when using this product.  
These may include, but are not limited to:
  - Dilated aortic root
  - Marked atrial enlargement
  - Inability to access the right atrium through the inferior vena cava
  - Rotated heart
  - Abnormal left atrial geometry
  - Vascular malformations

## Adverse Events

In addition to the complications associated with any cardiac catheterization, the following could occur during a transseptal catheterization:

- Inferior vena cava puncture
- Aorta puncture
- Atrial free wall puncture
- Coronary sinus puncture
- Arterial embolism from thrombus at the puncture site
- Tamponade
- Residual atrial septal defects
- Atrial arrhythmia

## Directions for Use

1. To determine the size and location of the left atrial septum, it may be helpful to perform right side angiography.
2. Thoroughly flush the transseptal needle.
3. Inspect all components prior to use for integrity and appropriateness for the particular procedure.
4. The transseptal needle/stylet set is advanced through the transseptal catheter and/or introducer until the tip of the needle is just within the catheter. Note: Ensure the transseptal needle is free to twist and/or rotate without resistance, as it is advanced to this position.
5. Retract the stylet prior to puncturing the interatrial septum.
6. The tip of the transseptal needle is positioned against the interatrial septum.
7. The transseptal needle is advanced through the septum into the left atrium once the correct location is confirmed using fluoroscopy (preferably biplane).
8. The successful puncture into the left atrium is confirmed by observing a left atrial pressure tracing. If an incorrect pressure tracing is observed, a small amount of contrast media is injected to identify the positioning.
9. The transseptal catheter is advanced over the transseptal needle into the left atrium.
10. The transseptal needle is slowly removed.

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	en
	Contents are non-pyrogenic and ethylene oxide sterilized.
	Do Not Resterilize
	Single use only. Do not reuse!
	Manufactured by
	Manufacturing Date
	Lot Number
	Use Before
	Contents
	Do not use if packaging is damaged!
	Catalog No./ Model No.
	Attention, see instructions
	Authorized European Representative
	Universal Product Number
	Proximal Gage
	Distal Gage
	This product is light sensitive. Do not use if stored outside protective outer carton. Store in a cool, dark, and dry place.

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Merit Medical Systems, Inc.
U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd,