Squirt™ FLUID DELIVERY SYSTEM

Contents of unopened, undamaged package are sterile and non-pyrogenic.
Contents: Squirt Fluid Delivery System

Caution: Federal (U.S.A) law restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE

Intended Use: The Squirt Syringe is intended for the controlled administration of thrombolytic agents into the peripheral vasculature.

Warning: This device must be used in accordance with appropriate labeling of thrombolytic solution and Fountain™ Infusion thrombolytic catheter or catheter of choice. See those package inserts for contraindications, side effects and precautions associated with thrombolytic therapy.

Using aseptic technique, peel open lid and remove Squirt from the tray.

Priming Instructions: Fill a syringe (not included) with sterile saline or with the medication of choice. (The Squirt is designed to accept either a 10 or 20 ml syringe.) Debubble the syringe using standard hospital protocol. This may include tapping the syringe with a hemostat or similar device.

Attach syringe to Squirt (see Figure 1). Make sure that the syringe connection is airtight. (The syringe rotator should be tightened by hand if using a syringe with a rotating adapter.)

Holding the Squirt in an upright position activate the trigger bar repeatedly until all air bubbles are out of the check valve area of the Squirt. (See Figure 1.)

This may include tapping the Squirt fluid path with a hemostat or similar device. [Note: Clinician should attach a small piece of tubing if concerned about fluid dripping out of the end of the Squirt during the priming process.]

Turn Squirt such that the Sherlock connector is pointing up. Activate the trigger bar until all air bubbles are out of fluid path. (See Figure 2). This may include tapping with a hemostat or similar device. This step may have to be repeated several times to fully debubble the system.

Filling with Thrombolytic Solution: If the clinician is desirous of conserving thrombolytic medication the Squirt should be primed as instructed above. After priming the system with saline the Squirt can be filled with a thrombolytic solution.


Turn the Squirt so the syringe is pointing down. (See Figure 3) Remove the priming reservoir syringe that is loaded with sterile saline. Fill a syringe with thrombolytic solution.

Using the thrombolytic solution syringe, place a small amount of thrombolytic solution into the female luer connector of the Squirt. This will cause a small meniscus of therapeutic solution to be placed on the female luer connector. (See Figure 4.)

Attach the syringe loaded with therapeutic solution. Make sure that the connection is airtight. If any bubbles inadvertently enter the system they may be removed by activating the trigger bar until all bubbles are out of the fluid path. (The dead space volume is approximately 0.5 ml.) The Squirt Fluid Delivery System is now ready to be attached to a catheter that is used for thrombolysis. **WARNING:** Always follow the thrombolysis catheter's instructions for use.

Adjusting the Stroke Volume of the Squirt: The stroke volume ejected from the Squirt can be adjusted from 0 - 1 ml of fluid by turning the knurled knob located in the handle. When holding the Squirt with the Sherlock connector pointed away from user, rotate the knob clockwise to decrease stroke volume. Rotating the knob in a counter-clockwise direction will increase the stroke volume of the device. (See figure 5.) Adjust the plunger tip to the amount of fluid to be infused with each stroke by aligning the plunger ring with the ml marks on the Squirt's barrel. After activating the trigger once the dosage is set and will deliver the same amount of fluid each time the activation trigger is fully pulled.

Instructions for use with a reservoir bag or bottle: Attach tubing connector to the female luer attachment located on the underside portion of the Squirt. Make sure the connection is airtight so that no air can enter the system. Holding the Squirt in an upright position prime the Squirt System in a similar manner to the priming instructions as listed above. Turn Squirt such that the Sherlock connector is pointing up. Repeatedly activate the trigger bar until all air is out of fluid path. This may include tapping the Squirt with a hemostat or similar device. The device is now ready to inject fluid into device or catheter of choice. (See figure 6) Caution: Make sure that fluid level in reservoir bag or bottle is continuously monitored so air is not injected inadvertently.

Patent Pending

MERIT MEDICA®

Manufacturer:
Merit Medical Systems, Inc.,
South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative:
Merit Medical Ireland Ltd
Parkmore Business Park West, Galway, Ireland

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