

Rad Board®

AND ACCESSORIES

Rad Trac™
RAD BOARD® SUPPORT BASE
MERTMEDICAL

Rad Board Xtra™
RAD BOARD® EXTENSION
MERTMEDICAL




MERTMEDICAL®

Efficiency and Protection

our radial advantage

The well-known benefits of the **RADIAL APPROACH** include less bleeding, quicker ambulation for patients, and lower procedure costs. Merit Medical offers a suite of radial-specific products designed to further optimize the **EFFICIENCY** of your radial procedures by streamlining and standardizing your set up. These radial accessories include a radial arm board with radiation scatter **PROTECTION** to enhance the **SAFETY** of clinical staff.¹



Rad Board®



The **Rad Board** is the only arm board that offers radiation scatter protection.

- 33-40% reduction at waist height¹
- 21-30% reduction at neck height¹

Provides quick, easy and standard set up for your radial access procedures

Fits all cath lab/radiology tables

Reversible and reusable

Uprights support drape and help form a reservoir for keeping devices and fluids on the board

Rad Rest®



Provides wrist and elbow support

Stabilizes and positions arm

Easily adjusts to accommodate any upper extremity length

Soft foam for increased patient comfort

Rad Board Xtra™



Allows 90° perpendicular extension of arm during access

Provides an additional 11" (28 cm) of working space when placed next to Rad Board

Reversible and reusable

THINK RADIAL. **DISCOVER MERIT™**

From set up to hemostasis, Merit Medical offers a complete radial product portfolio to help you achieve the procedural **EFFICIENCY** and **SAFETY** that is essential for delivering the best possible **CARE**.

ORDERING INFORMATION

Part #	Description	Order Unit
RB100	Rad Board® arm board	Sold individually, non-sterile
RBX	Rad Board Xtra™	Sold individually, non-sterile
RT913	Rad Trac™	Sold individually, non-sterile
RR811	Rad Rest®	Individually packaged, 10 units per box, non-sterile

¹Data on file.

Rad Trac™



Allows easy placement of the Rad Board when patient is on the table

Fits all cath lab/radiology tables

Reusable



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
www.merit.com



Merit Medical Europe, Middle East, and Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country

Austria 0800 295 374 | Belgium 0800 72 906 (Dutch)
0800 73 172 (French) | Denmark 80 88 00 24
France 0800 91 60 30 | Finland 0800 770 586 |
Germany 0800 182 0871 | Ireland (Republic) 1800 553 163
Italy 800 897 005 | Luxembourg 8002 25 22
Netherlands 0800 022 81 84 | Norway 800 11629
Sweden 020 792 445 | UK 0800 973 115