



## *The Right Combination*

*Everything necessary for the procedure is on  
your shelf in an organized, ready-to-use tray*

PERICARDIOCENTESIS KIT™

# BE PREPARED!

## THE RIGHT COMBINATION

Everything necessary for the procedure is on your shelf in an organized, ready-to-use tray.

## TIME SAVINGS

Time previously spent assembling a kit from scratch can now be focused toward the patient.

## COST SAVINGS

Ordering the ready-to-use kit saves on personnel time and lowers overall component cost.

## PROCEDURE SPECIFIC DESIGN

High-flow catheters are developed with pericardiocentesis procedures in mind.



### 8.3F KIT COMPONENTS:

- 8.3F (40 cm) pigtail or straight catheter (with marker band)
- 8.5F (22 cm) dilator
- .035" (0.89 mm) x 80 cm, double ended (J & straight) standard PTFE guide wire
- 25G x 1 1/2" (3.8 cm) anesthesia needle
- 18G x 1 1/2" (3.8 cm) medication needle
- 18G x 3.5" (9 cm) access needle
- 18G x 5.9" (15 cm) access needle
- #11 Futura® safety scalpel
- 10 mL luer lock syringe
- 20 mL luer lock syringe
- 60 mL luer lock syringe
- 1000 mL disposal bag with drainage spout
- (2) Large bore three-way stopcocks
- Double-ended alligator clip cable
- (10) 4" x 4" (10 cm x 10 cm) gauze sponges
- 60 cm x 60 cm fenestrated drape
- 46 cm x 60 cm absorbent barrier tri-layer drape
- Suture 3-0 silk with curved cutting needle
- PAL™ label set – (4) yellow saline, (4) red specimen, (4) blue lidocaine
- 2 ChloroPrep®\* 3 mL (orange tint) for use in PC801 and PC802



### 6F KIT COMPONENTS:

- 6F (60 cm) high-flow pigtail or straight catheter
- 6F (22 cm) dilator
- .035" (0.89 mm) x 80 cm, double ended (J and straight) heavy duty PTFE guidewire
- 22G x 2 3/8" (6 cm) anesthesia needle
- 19G x 1 1/2" (3.8 cm) medication needle
- 18G x 3.5" (9 cm) access needle
- 18G x 5.9" (15 cm) access needle
- #11 scalpel
- 10 mL luer lock syringe
- 60 mL luer lock syringe
- 1400 mL disposal bag with connecting tube
- Three-way stopcock
- Double-ended alligator clip cable
- (10) 4" x 4" (10 cm x 10 cm) gauze sponges
- 40 cm x 40 cm fenestrated drape
- 46 cm x 60 cm absorbent barrier tri-layer drape

\* ChloroPrep® is a registered trademark of CareFusion 2200, Inc.

Cat. No.	Description
----------	-------------

PC801	8.3F Pigtail Catheter
PC802	8.3F Straight Catheter
PC803EU	8.3F Pigtail Catheter, Europe
PC804EU	8.3F Straight Catheter, Europe
PC101	6F Pigtail Catheter
PC201	6F Straight Catheter

