



# KanguruWeb™

*Every patient deserves your best care*

# KanguruWeb™

Every patient deserves your best care

The KanguruWeb system provides consistent, reliable retraction of the abdominal pannus. The disposable contoured web eliminates potential skin complications that can occur with using tape. Use the lockable arm rests on all patients and the disposable web as needed for obese patients.

## ORDERING INFORMATION

Catalog numbers:		Packaged
<b>DW-3600</b>	Disposable Web	25 per box
<b>KANG100</b>	Kanguru Access Kit, Sterile	10 per box
	(1) Disposable web	
	(1) 18G x 9cm Needle	
	(1) 22G x 6cm needle	
<b>RB-2500</b>	Fluoroscopy Base – Reusable	1 per box
<b>LA-8600</b>	Lockable Arm Rest – Reusable	2 per box

## FEATURES INCLUDE:

**Abdominal Pannus Relocation** - Improved femoral site access

**Rapid Device Deployment** - Easier and expedient patient preparation

**Preserves Dignity and Respect** - Promotes a feeling of security and comfort for patients


**Facilitates Skin Integrity** - Helps eliminate complications from using tape - no skin burns, tears or abrasions

**Disposable Web Design** - Quickly and easily retracts the abdominal pannus in less than one minute


**Flexibility** - Easily adjusts to accommodate different body shapes without changing the patient prep time, regardless of pannus shape or size

**Lockable Arm Rests** - Prevents arms from falling off the table / decreases potential contamination


**Radiolucent Base** - Prevents fluoroscopy / imaging interference



**PRE WEB PLACEMENT**  
Lateral view, hidden femoral access site.



**PRE WEB DEPLOYMENT**  
Frontal view, hidden femoral access site.



**POST WEB DEPLOYMENT**  
Lateral view, exposed femoral access site.

**PACKAGED FOR YOUR CONVENIENCE**



### KANGURU ACCESS KIT

Disposable web  
18G x 9cm puncture needle  
22G x 6cm anesthesia needle

