



MERITMEDICAL | Inflation Devices & Accessories

basixTOUCH™

ENHANCED ANALOG INFLATION DEVICE

Fast and versatile – the basixTOUCH

provides expanded capacity for inflation. With 35 atm pressure capability and 30 mL syringe volume, the basixTOUCH is ideal for any interventional lab using high pressure, large volume balloons. With wider threads and quick-release handle, the basixTOUCH offers rapid inflation and deflation to busy labs that need top efficiency.



Catalog #	Product Description (5 per box)
IN8130	BasixTOUCH
IN8112	BasixTOUCH with Access-9
IN8152	BasixTOUCH with AccessPLUS
IN8302	BasixTOUCH with Honor
IN8352	BasixTOUCH with Honor + Extension Tubing
IN8403	BasixTOUCH with MBA

BLUEDIAMOND™

DIGITAL INFLATION DEVICE + FLUID DISPENSING SYRINGE

With a precise digital pressure gauge

and digital timer, you don't need to wonder how long an inflation was, or what pressure was obtained. The Blue Diamond digital inflation device offers an integrated digital gauge and timer, with the ability to record and recall prior inflation information. Ideal for ensuring precise inflation pressure and for tracking procedure data.

Catalog #	Product Description (5 per box)
IN7130	Blue Diamond
IN7112	Blue Diamond with Access-9
IN7152	Blue Diamond with AccessPLUS
IN7352	Blue Diamond with Honor + Extension Tubing
IN7403	Blue Diamond with MBA



For 25 years, Merit Medical has been the market leader in inflation device technology. By working closely with clinicians, we are able to create innovative devices to meet the need for faster inflation, greater capacity, accuracy, and reliability. Our broad product portfolio offers a wide range of features and price options—analogue and digital, low and high pressure capacity—for inflation, deflation, and measurement of pressures in a variety of balloons and stents.

basixCOMPAK™

ANALOG INFLATION DEVICE

The gold standard inflation device for

interventional balloon inflation – with over 10 million devices used worldwide, busy interventional labs can rely on the basixCOMPAK for their procedures. The basixCOMPAK offers 30 atm pressure and 20 mL volume and can be packaged with a wide variety of angioplasty kits and procedure packs.



Catalog #	Product Description (5 per box)
IN4130	BasixCOMPAK
IN4230	BasixCOMPAK with Access-9
IN4530	BasixCOMPAK with AccessPLUS
IN4330	BasixCOMPAK with Honor
IN4352	BasixCOMPAK with Honor + Extension Tubing
IN4430	BasixCOMPAK with MBA

ENDOTEK BIG60®

LARGE VOLUME ANALOG INFLATION DEVICE

The BIG60 inflation device is designed

to rapidly inflate and deflate large volume, low pressure balloons. With 60 mL fluid volume and a 12 atm pressure capability, the BIG60 is ideal for interventional practitioners that routinely use large volume interventional devices.



Catalog #	Product Description (5 per box)
ENDO-AN6012	BIG60



HEMOSTASIS VALVES AND ACCESSORIES

Patients and medical staff want to reduce blood loss as much as possible.

In response, Merit Medical has developed a full line of hemostasis valves to minimize blood loss in interventional medicine, while supporting a wide array of procedures and devices. Merit Medical has been a market leader in the design and manufacture of hemostasis valves and has manufactured over 10 million valves to support the diverse needs of our users.



Access-9™

Premier large bore 9F device with standard Tuohy seal with fine pitch threads for rapid opening and closure.

Catalog # Access-9 Accessories

MAP111	Metal insertion tool
MAP112	Metal insertion tool, torque device



MAP500
Torque Device



AccessPLUS™

Premier large bore 9F device with specially designed Tuohy valve that improves wire grip.

Catalog # AccessPLUS Accessories

MAP150	Standalone
MAP151	Metal insertion tool
MAP152	Metal insertion tool, torque device



MAP550
Metal Insertion Tool



DoublePlay™

Large bore double valve for kissing balloon procedures

Catalog # DoublePlay Accessories

MAP220	Metal insertion tool
MAP222	Metal insertion tool, torque device (2)



MAP600
Plastic Insertion Tool



Honor®

Small bore 7F device with push-pull bleedback mechanism. The Honor has an innovative triple-slit bleedback seal, and no Tuohy-style seal.

Catalog # Honor Accessories

MAP302	Metal insertion tool, torque device
MAP352	Plastic insertion tool, 20 cm bonded extension tubing, torque device
MAP353	Metal insertion tool, 50 cm bonded extension tubing, torque device



MBA™

Small bore device combines both the reliable Tuohy seal as well as a bleedback seal, providing operators with 2 hemostatic mechanisms.

Catalog # MBA Accessories

MAP400	Standalone
MAP403	Metal insertion tool, torque device
MAP453	Metal insertion tool, 10 cm bonded extension tubing, torque device



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
www.merit.com

Merit Medical Europe, Middle East, and Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country

Austria 0800 295 374 | Belgium 0800 72 906 (Dutch)
0800 73 172 (French) | Denmark 80 88 00 24
France 0800 91 60 30 | Finland 0800 770 586 |
Germany 0800 182 0871 | Ireland (Republic) 1800 553 163
Italy 800 897 005 | Luxembourg 8002 25 22
Netherlands 0800 022 81 84 | Norway 800 11629
Sweden 020 792 445 | Switzerland 800 9000 5060
UK 0800 973 115

ID 092214 400899001/B