

CONCIERGE™

— GUIDING CATHETER —

Quality and Performance...At Your Service



Merit's advanced braiding technology and proprietary polymer blend shaft materials combined with progressively compliant distal segments allow for an increase in lumen size while maintaining the exceptional **back-up support** afforded by the Concierge product line.

The catheter tip responds directly to **torque** applied to the proximal section of the catheter.


Optimal stiffness at room temperature allows for **better pushability** and makes the catheter less prone to kinking—providing superior insertability. The catheter incorporates a thermally responsive proprietary polymer blend that provides **better compliance** at body temperature.


Incorporating an atraumatic **soft tip**, the catheter reaches the coronary ostium with **accuracy** and ease.



Concierge Catalog Numbers by French Size

Shape	Sideholes	Tip	5F	6F	7F	8F
Judkins Right	-	3.0	CGC5100JR30	CGC6100JR30	CGC7100JR30	CGC8100JR30
Judkins Right	-	3.5	CGC5100JR35	CGC6100JR35	CGC7100JR35	CGC8100JR35
Judkins Right	2	3.5	CGC5100JR35SH	CGC6100JR35SH	CGC7100JR35SH	CGC8100JR35SH
Judkins Right	-	4.0	CGC5100JR40	CGC6100JR40	CGC7100JR40	CGC8100JR40
Judkins Right	2	4.0	CGC5100JR40SH	CGC6100JR40SH	CGC7100JR40SH	CGC8100JR40SH
Judkins Right	-	5.0	CGC5100JR50	CGC6100JR50	CGC7100JR50	CGC8100JR50
Judkins Right	2	5.0	CGC5100JR50SH	CGC6100JR50SH	CGC7100JR50SH	CGC8100JR50SH
Judkins Right	-	6.0	CGC5100JR60	CGC6100JR60	CGC7100JR60	CGC8100JR60
Judkins Left	-	3.0	CGC5100JL30	CGC6100JL30	CGC7100JL30	CGC8100JL30
Judkins Left	-	3.5	CGC5100JL35	CGC6100JL35	CGC7100JL35	CGC8100JL35
Judkins Left	2	3.5	CGC5100JL35SH	CGC6100JL35SH	CGC7100JL35SH	CGC8100JL35SH
Judkins Left (ST)	-	3.5	CGC5100JL35ST	CGC6100JL35ST	CGC7100JL35ST	CGC8100JL35ST
Judkins Left	-	4.0	CGC5100JL40	CGC6100JL40	CGC7100JL40	CGC8100JL40
Judkins Left	2	4.0	CGC5100JL40SH	CGC6100JL40SH	CGC7100JL40SH	CGC8100JL40SH
Judkins Left (ST)	-	4.0	CGC5100JL40ST	CGC6100JL40ST	CGC7100JL40ST	CGC8100JL40ST
Judkins Left	-	4.5	CGC5100JL45	CGC6100JL45	CGC7100JL45	CGC8100JL45
Judkins Left	2	4.5	CGC5100JL45SH	CGC6100JL45SH	CGC7100JL45SH	CGC8100JL45SH
Judkins Left	-	5.0	CGC5100JL50	CGC6100JL50	CGC7100JL50	CGC8100JL50
Judkins Left	2	5.0	CGC5100JL50SH	CGC6100JL50SH	CGC7100JL50SH	CGC8100JL50SH
SBS 3.0	-	3.0	CGC5100SBS30	CGC6100SBS30	CGC7100SBS30	CGC8100SBS30
SBS 3.0	2	3.0	CGC5100SBS30SH	CGC6100SBS30SH	CGC7100SBS30SH	CGC8100SBS30SH
SBS 3.5	-	3.5	CGC5100SBS35	CGC6100SBS35	CGC7100SBS35	CGC8100SBS35
SBS 3.5	2	3.5	CGC5100SBS35SH	CGC6100SBS35SH	CGC7100SBS35SH	CGC8100SBS35SH
SBS 3.5 (Mod)	-	3.5	CGC5100SBS35-M	CGC6100SBS35-M	CGC7100SBS35-M	CGC8100SBS35-M
SBS 4.0	-	4.0	CGC5100SBS40	CGC6100SBS40	CGC7100SBS40	CGC8100SBS40
SBS 4.0	2	4.0	CGC5100SBS40SH	CGC6100SBS40SH	CGC7100SBS40SH	CGC8100SBS40SH
SBS 4.0 (Mod)	-	4.0	CGC5100SBS40-M	CGC6100SBS40-M	CGC7100SBS40-M	CGC8100SBS40-M






P R O D U C T D E S C R I P T I O N			
5F	6F	7F	8F
G U I D E O D			
0.070"(1.8mm)	0.082"(2.1mm)	0.092"(2.3mm)	0.105"(2.7mm)
G U I D E I D			
0.057"(1.4mm)	0.070"(1.8mm)	0.078"(2.0mm)	0.088"(2.2mm)
Braiding- Stainless Steel Inner Lining - PTFE (polytetrafluoroethylene) Outer Coating - (proprietary polymer blend)			

Concierge Catalog Numbers by French Size

Shape	Sideholes	Tip	5F	6F	7F	8F
Amplatz Right 1	-	AR1	CGC5100AR1	CGC6100AR1	CGC7100AR1	CGC8100AR1
Amplatz Right 1	2	AR1	CGC5100AR1SH	CGC6100AR1SH	CGC7100AR1SH	CGC8100AR1SH
Amplatz Right 2	-	AR2	CGC5100AR2	CGC6100AR2	CGC7100AR2	CGC8100AR2
Amplatz Right 2	2	AR2	CGC5100AR2SH	CGC6100AR2SH	CGC7100AR2SH	CGC8100AR2SH
Amplatz Right 3	-	AR3	CGC5100AR3	CGC6100AR3	CGC7100AR3	CGC8100AR3
Amplatz Left 1	-	AL1	CGC5100AL1	CGC6100AL1	CGC7100AL1	CGC8100AL1
Amplatz Left 1	2	AL1	CGC5100AL1SH	CGC6100AL1SH	CGC7100AL1SH	CGC8100AL1SH
Amplatz Left 2	-	AL2	CGC5100AL2	CGC6100AL2	CGC7100AL2	CGC8100AL2
Amplatz Left 2	2	AL2	CGC5100AL2SH	CGC6100AL2SH	CGC7100AL2SH	CGC8100AL2SH
Amplatz Left 3	-	AL3	CGC5100AL3	CGC6100AL3	CGC7100AL3	CGC8100AL3
Multipurpose 1	-	MP	CGC5100MP1	CGC6100MP1	CGC7100MP1	CGC8100MP1
Multipurpose 1	2	MP	CGC5100MP1SH	CGC6100MP1SH	CGC7100MP1SH	CGC8100MP1SH
Left Coronary Bypass	-	LCB	CGC5100LCB	CGC6100LCB	CGC7100LCB	CGC8100LCB
Left Coronary Bypass	2	LCB	CGC5100LCBSH	CGC6100LCBSH	CGC7100LCBSH	CGC8100LCBSH
Right Coronary Bypass	-	RCB	CGC5100RCB	CGC6100RCB	CGC7100RCB	CGC8100RCB
Right Coronary Bypass	2	RCB	CGC5100RCBSH	CGC6100RCBSH	CGC7100RCBSH	CGC8100RCBSH
Internal Mammary	-	IM	CGC5100IM	CGC6100IM	CGC7100IM	CGC8100IM
Internal Mammary	2	IM	CGC5100IMSH	CGC6100IMSH	CGC7100IMSH	CGC8100IMSH
El Gamal	-	EG	CGC5100EG	CGC6100EG	CGC7100EG	CGC8100EG
Hockey Stick	-	HS	CGC5100HS	CGC6100HS	CGC7100HS	CGC8100HS



All catheters available in 100cm length.

Other Products from Merit Medical



Merit's market-leading analog and digital inflation devices provide reliable and accurate inflation for your most demanding cases. Ergonomically designed with accuracy and quality in mind, Merit inflation devices are also available in custom kit configurations.



Merit Medical is pleased to offer a variety of safety and waste related products. These products have been designed to complement other devices in an effort to make the clinical setting safer for both the clinician and the patient.



Merit has developed a complete line of technically sophisticated, clinically acclaimed hemostasis valves and angioplasty accessories.



Merit diagnostic cardiology catheters are available in Softouch and Performa styles. Designed with a nylon shaft and wire braid, Softouch catheters are the original workhorse catheters of Merit Medical. The Performa incorporates an improved nylon shaft and flatwire braid.



Merit is proud to offer an increased level of service for customers through direct tray and pack distribution. To streamline this new service offering, Merit has designed a variety of trays and packs which are readily available to meet various clinical needs.



Merit Medical offers a broad line of introducers and accessories. Product offerings include the Prelude® Sheath Introducers and accessories, Merit MAK™ and S-MAK™ Introducers and accessories.



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, AFRICA AND THE MIDDLE EAST (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone for specific country:

Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 Germany 0800 182 08 71

Ireland (republic) 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

400912001/B ID 010610