



OUTAKE™

*The Merit Tunneled Catheter
Extraction Kit is Convenient,
Simple, and Cost Effective.*



OUTAKE™

The OuTake Catheter Extractor is a convenient, simple, and cost effective device used to remove tunneled double lumen hemodialysis catheters with an ingrown cuff.

CONVENIENT

The OuTake Kit contains all necessary components to perform a catheter extraction procedure at the bedside, out-patient procedure room, or in the office.

There is no need to gather components from other departments or resterilize. The procedure can be performed by a Physician Assistant or Nurse Practitioner.

SIMPLE

The OuTake Kit promotes a streamlined procedure with less trauma to the patient. The OuTake device enters the exit hole and cuts ingrown tissue at the same time to facilitate smooth catheter removal.

WHY USE THE OUTAKE?

- *Reduces bleeding compared to traditional procedures*
- *Reduces likelihood of scar formation*
- *Potentially reduces the need for follow up*
- *Decreases procedure time and cost*
- *Improves patient comfort*

The Catheter Extractor is indicated for use with tunneled, double lumen hemodialysis catheters with cuff sizes ranging from 10.8 French to 14.5 French.

COST EFFECTIVE

The OuTake Kit reduces the overall cost of the procedure by having all components in one location, eliminating lab time, and shortening the length of the procedure. Ordering with one catalog number reduces inventory and increases shelf space.

*OuTake™ device
not to scale with other
items in this kit*



KIT CONTENTS



CET100 - Stand Alone OuTake™

CET100 - Stand Alone OuTake™ Catheter Extractor

CET200 - Catheter Extraction Kit with OuTake™

- (1) OuTake™ Catheter Extractor
- (1) Lidocaine, 2%
- (1) ND. 18G x 1 ½
- (1) ND. 25G x 1 ½
- (1) 24 x 24 Poly Wrap
- (1) Oint. Triple Antibiotic, 1/32oz.
- (1) Shortstop Sharps Holder
- (1) Safety Scalpel, #11
- (1) ChloroPrep, 10.5 ml
- (1) Tegaderm, 10x12cm Framed
- (3) Sterile Strips, ½ x 4"
- (1) Syringe, 12cc Luer Lock
- (1) Drape, 18x26 Oval Fenestration W/Adhesive
- (1) Iris Scissors, Straight, 4 ½"
- (1) Mosquito, Curved 5"
- (10) Gauze 4x4, 12ply
- (2) Cloth Towels, Blue
- (1) Gloves, Nitrile Exam, Medium
- (1) Transport Tube, 10ml Clear W/ Screw Cap



CET200 - Catheter Extraction Kit with OuTake™

CET300 - Catheter Extraction Kit without OuTake™

- (2) Lidocaine, 2%
- (2) ND. 18G x 1 ½
- (2) ND. 25G x 1 ½
- (2) 24 x 24 Poly Wrap
- (2) Oint. Triple Antibiotic, 1/32oz.
- (2) Shortstop Sharps Holder
- (2) Safety Scalpel, #11
- (2) ChloroPrep, 10.5 ml
- (2) Tegaderm, 10x12cm Framed
- (3) Sterile Strips, ½ x 4"
- (1) Syringe, 12cc Luer Lock
- (1) Drape, 18x26 Oval Fenestration W/Adhesive
- (1) Iris Scissors, Straight, 4 ½"
- (1) Mosquito, Curved 5"
- (10) Gauze 4x4, 12ply
- (2) Cloth Towels, Blue
- (1) Gloves, Nitrile Exam, Medium
- (1) Transport Tube, 10ml Clear W/ Screw Cap



CET300 - Catheter Extraction Kit without OuTake™

