



« 10F »
NOW AVAILABLE

EXCEPTIONAL INFLATION CONTROL



Q50[®]

Stent Graft Balloon Catheters

The aortic modeling balloons that offer the broadest diameter range with the lowest available profile.

Q50

Stent Graft Balloon Catheters

OPTIONS WITH PERFORMANCE

Q50 stent graft balloon catheters have demonstrated reliability for abdominal and thoracic EVAR procedures and reinterventions. In addition to reliable stent graft modeling and expansion, the Q50 balloon catheters are ideal for temporary occlusion of large vessels.

PRODUCT FEATURES

- Broadest balloon diameter range with the lowest available profile
 - 10mm to 50mm inflation diameters
 - 10F and 12F introducer sheath compatible
- Strong, yet compliant balloon material
- Two versatile catheter lengths, 65 cm and 100 cm
- Unique, tapered flexible tip facilitates device passage
- Exceptional balloon inflation and deflation rates
- Optimized catheter pushability


BALLOON COMPLIANCE CHART

Recommended Inflation Volume (cc)	Balloon Diameter (mm)
3	10
6	20
16	30
32	40
44	45
60	50

ORDERING INFORMATION

Part Number	Sheath Compatibility	Balloon Diameter Range	Catheter Length
Q50X™ ●			
Q50-65-X	10F	10–50 mm	65 cm
Q50-100-X	10F	10–50 mm	100 cm
Q50® PLUS ●			
Q50-65P	12F	10–50 mm	65 cm
Q50-100P	12F	10–50 mm	100 cm

Manufactured by
QXMedical



Q50; Q50X™ and Q50® PLUS are trademarks of QXMedical, LLC

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115