




True Form™

RESHAPABLE GUIDE WIRE


- Shapeable tip retains shape during procedures
- Shapeable tip may be shaped and reshaped multiple times for vessel cannulation¹
- Stainless Steel core gives excellent support for devices
- Hydrophilic coating to increase trackability through vessels
- Flexible shaft to navigate tortuous anatomy



SHAPEABLE TIP TECHNOLOGY THAT TRULY HOLDS ITS SHAPE


True Form™

RESHAPABLE GUIDE WIRE



GUIDE WIRE DRAG

(Frictional Force (g))



The True Form hydrophilic coating is up to 7 times more lubricious than the competition.¹
 Minimum of 15 samples tested

“The ability of the True Form Guide Wire to be shaped multiple times and keep its shape allows me to use the same guide wire, decreasing the need to use additional guide wires.”¹



Dr. Shinichi Hori MD, PhD.
 Chief Physician and Director,
 Gate Tower Institute for Image Guided Therapy
 Dr. Hori is a paid consultant of Merit Medical Systems, Inc.

Product Number	OD	Length (cm)	Tip Shape	Distal Taper (cm)	Shapeable Tip (cm)	Distal Coil Length (cm)	Tip Load (g)
TF14145S	0.014"	145	Straight	36	2	5	3.3
TF14165S	0.014"	165	Straight	36	2	5	3.3
TF14180S	0.014"	180	Straight	36	2	5	3.3
TF14180A	0.014"	180	Angled	36	2	5	3.3

1. Data on file

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
 1600 West Merit Parkway
 South Jordan, Utah 84095
 1.801.253.1600
 1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
 Amerikalaan 42, 6199 AE Maastricht-Airport
 The Netherlands
 +31 43 358 82 22

Merit Medical Ireland Ltd.
 Parkmore Business Park West
 Galway, Ireland
 +353 (0) 91 703 733

merit.com

Australia
 +61 (0) 3 8373 4943

Austria
 0800 295 374

Belgium
 0800 72 906 (Dutch)
 0800 73 172 (Français)

Brazil
 +55 11 4561 - 0062

Canada
 1.800.364.4370

China
 +86-10 8561 0788

Denmark
 80 88 00 24

Finland
 0800 770 586

France
 0800 91 60 30

Germany
 0800 182 0871

Hong Kong
 +852-2207 6888

India
 +91-80-41223376

Ireland (Republic)
 1800 553 163

Italy
 800 897 005

Luxembourg
 8002 25 22

Netherlands
 0800 022 81 84

Norway
 800 11629

Portugal
 308 801 034

Russia
 +7 495 221 89 02

Singapore
 +65 67536311

South Korea
 +8210 2205 0707

Spain
 +34 911238406

Sweden
 020 792 445

Switzerland
 (Deutsch)
 +41 22 518 02 30

(Français)
 +41 22 518 02 52

(Italiano)
 +41 22 518 00 35

UK
 0800 973 115