


RADIAL APPROACH

IT'S TIME FOR OPTIONS



*Begin procedures
with a Prelude®*

Prelude[®] to Finale[™]

Prelude[®]

Merit Advance[™]
ANGIOGRAPHIC NEEDLES

InQwire[®]

CONCIERGE[™]
GUIDING CATHETER

Performa[®]

SHEATH INTRODUCERS

Excellent transitions facilitate smooth entry into patient and vessel. Resists kinking during insertion.



ANGIOGRAPHY NEEDLES

Ultra sharp needle bevel and lubricious cannula ensures smooth, easy, atraumatic vascular access. Ergonomic hub and bevel indicator enhance clinical confidence.



GUIDE WIRES

InQwire[®] product line represents a complete offering of standard, fixed and moveable core guide wires. Hydrophilic guide wires are also available in a wide range of configurations.



GUIDE CATHETERS

Advanced braiding technology and proprietary shaft material provide the clinician exceptional torque and back-up support. Optimal stiffness allows for better pushability while maintaining an atraumatic soft tip. (Not yet available in U.S.)



DIAGNOSTIC CARDIOLOGY CATHETERS

The Ultimate[™] line of diagnostic angiography catheters for procedures performed using the radial approach.



Merit is pleased to offer several new product options in a comprehensive radial approach product line. Combining Merit's broad offering with your clinical expertise can help achieve the improved outcomes you expect from the radial approach:

Reduced Complications

Quicker Ambulation

Increased Patient Comfort

Merit . . . your partner for Radial Procedures



INJECT 8[®]



MERIT MISER[®]

Finale[®]

HEMOSTASIS VALVES

Multiple options and sizes to meet your procedural demands and preferences.



SYRINGES

Generate more pressure with less applied force through smaller catheters using Inject8 and Inject10n[™] coronary control syringes.



INFLATION DEVICES

The market leader is proud to offer information at the touch of a button. A portable digital inflation device for the clinician who demands accuracy and convenience.



FLUID MANAGEMENT

Miser[®] contrast management systems, Meritrans[®] physiological pressure transducers, manifolds, stopcocks, and high pressure tubing are just a few of our comprehensive line of fluid management products.



RADIAL COMPRESSION

The Finale[™] and RadStat[®] radial compression devices provide controlled hemostasis, comfortable immobilization (RadStat) and more rapid ambulation after radial artery puncture.



End procedures with a Finale[™]



Founded in 1987, Merit Medical Systems is a world leader in the development, manufacture and distribution of proprietary disposable medical products used in interventional and diagnostic procedures, particularly in cardiology and radiology.

The company serves client hospitals worldwide with a domestic and international sales force totaling over 125 individuals.

For ordering information call Merit Medical customer service 1-800-35-MERIT.



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40

Free phone for specific country:

Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71

Republic of Ireland 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

401911001/A ID 083110