

ReSolve[®]

BILIARY LOCKING DRAINAGE CATHETER



MERITMEDICAL[®] | Targeted Drainage Solutions

ReSolve[®]

BILIARY LOCKING DRAINAGE CATHETER

Clinical Benefits

Polyurethane material allows for **long-term indwelling** time

Large drainage holes help **maximize drainage**


Catheter is designed to allow **unimpeded** drainage

Depth markers and radiopaque marker band helps with **fast and accurate** positioning

Atraumatic pigtail designed to reduce the potential for damaging the duodenum

Pigtail straightener allows for **easy insertion** of the stiffening cannulas


Hydrophilic coating facilitates **smooth insertion**





*Two drainage hole patterns
to accommodate multiple
lesion locations and varying
patient anatomy.*



Merit has up to 4X more drainage area



ReSolve[®] BILIARY LOCKING DRAINAGE CATHETER



RBC Ordering Information

RBDC Ordering Information

Catalog Number	French Size	Drainage Location	Pigtail to Distal Shaft Drainage Hole	Drainage Holes
RBC-8-038	8.5	14 cm	5 cm	17
RBC-10-038	10	14 cm	5 cm	17
RBC-12-038	12	14 cm	5 cm	18
RBC-14-038	14	14 cm	5 cm	18

Catalog Number	French Size	Drainage Location	Pigtail to Distal Shaft Drainage Hole	Drainage Holes
RBDC-8-038	8.5	11 cm	2 cm	17
RBDC-10-038	10	11 cm	2 cm	17
RBDC-12-038	12	11 cm	2 cm	18
RBDC-14-038	14	11 cm	2 cm	18

with Revolution[™] Catheter Securement Device

RBC-8-REV	8.5	14 cm	5 cm	17
RBC-10-REV	10	14 cm	5 cm	17
RBC-12-REV	12	14 cm	5 cm	18
RBC-14-REV	14	14 cm	5 cm	18

with Revolution Catheter Securement Device

RBDC-8-REV	8.5	11 cm	2 cm	17
RBDC-10-REV	10	11 cm	2 cm	17
RBDC-12-REV	12	11 cm	2 cm	18
RBDC-14-REV	14	11 cm	2 cm	18

with StayFIX[®] Fixation Device

RBC-8-SFX	8.5	14 cm	5 cm	17
RBC-10-SFX	10	14 cm	5 cm	17
RBC-12-SFX	12	14 cm	5 cm	18
RBC-14-SFX	14	14 cm	5 cm	18

with StayFIX Fixation Device

RBDC-8-SFX	8.5	11 cm	2 cm	17
RBDC-10-SFX	10	11 cm	2 cm	17
RBDC-12-SFX	12	11 cm	2 cm	18
RBDC-14-SFX	14	11 cm	2 cm	18

Packaged 5 per box

Contact Merit Medical for more information on our line of drainage products.



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
www.merit.com

Merit Medical Europe, Middle East,
and Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country

Austria 0800 295 374 | Belgium 0800 72 906 (Dutch)
0800 73 172 (French) | Denmark 80 88 00 24
France 0800 91 60 30 | Finland 0800 770 586 |
Germany 0800 182 0871 | Ireland (Republic) 1800 553 163
Italy 800 897 005 | Luxembourg 8002 25 22
Netherlands 0800 022 81 84 | Norway 800 11629
Sweden 020 792 445 | Switzerland 800 9000 5060
UK 0800 973 115

ID 031714 402324001/B