

# BiliQUICK™

## CHOLANGIOGRAPHY INJECTION KIT



 **MERTMEDICAL**  
ENDOTEK®

# BiliQUICK™ Cholangiography Rapid Refill Continuous Injection Kit



- Attach contrast, prime tubing, and connect device with the **All-In-One Design**
- Long, thin barrel **reduces required push force**
- **CRV Valve** allows pushing and pulling in a continuous motion without needing to refill individual syringes
- **Large Bore Tubing** allows smooth passage of thick contrast fluids

Catalog Number	Description	Packaged
K12-05301	Cholangiography Injection Kit	5 units



Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle  
East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733