

Embosphere®


Microspheres

EmboGold®

Microspheres

CATHETER COMPATIBILITY CHART


Merit Maestro®
Microcatheter



Merit Pursue™
Microcatheter



SwiftNINJA®
Steerable Microcatheter



		50-100	40-120	100-300	300-500	500-700	700-900	900-1200
Catheter	Inner Diameter							
Catheter	0.038" (0.97 mm)							
Catheter	0.035" (0.89 mm)							
EmboCath® Plus Microcatheter	0.028" (0.71 mm)							
Merit Maestro® Microcatheter 2.9F	0.027" (0.69 mm)							
Merit Maestro Microcatheter 2.8F	0.024" (0.62 mm)							
SwiftNINJA® Steerable Microcatheter 2.4F	0.021" (0.54 mm)							
Merit Maestro Microcatheter 2.4F	0.020" (0.52 mm)							
Merit Pursue™ Microcatheter 2.0F	0.020" (0.50 mm)							
Merit Maestro Microcatheter 2.1F	0.018" (0.46 mm)							
Merit Pursue Microcatheter 1.7F	0.016" (0.40 mm)							

Above catheter compatibility information based on internal testing.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733



Biosphere Medical, S.A.
Parc des Nations - Paris Nord 2
383 rue de la Belle Etoile
95700 Roissy en France
France

