

Embosphere®

Microspheres

EmboGold®

Microspheres

CATHETER COMPATIBILITY CHART

The elastic properties of Embosphere and EmboGold microspheres allow them to temporarily compress up to 33%, allowing for smooth microcatheter passage. Once through the microcatheter, the microspheres return to their original spherical shape and stated diameter for predictable, targeted delivery.

Compatible Microsphere Size (µm)

Catheter	Inner Diameter	50-100	40-120	100-300	300-500	500-700	700-900	900-1200
Catheter	0.038" (0.97 mm)							
Catheter	0.035" (0.89 mm)							
EmboCath® Plus Microcatheter	0.028" (0.71 mm)							
Merit Maestro® Microcatheter	0.027" (0.69 mm)							
Merit Maestro Microcatheter	0.024" (0.61 mm)							
Microcatheter	0.023" (0.58 mm)							
SwiftNINJA® Steerable Microcatheter	0.021" (0.54 mm)							
Merit Maestro Microcatheter	0.020" (0.51 mm)							
Microcatheter	0.019" (0.48 mm)							
Merit Maestro Microcatheter	0.018" (0.46 mm)							
Microcatheter	0.017" (0.43 mm)							
Microcatheter	0.016" (0.41 mm)							

Microspheres depicted are not to scale and are enlarged for clarity. Above catheter compatibility information based on internal testing.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30
(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115