

DESIGNED FOR COMFORT & ACCURACY



ReSolve®
LOCKING DRAINAGE CATHETER

ReSolve®

LOCKING DRAINAGE CATHETER

To prevent irritation of the surrounding tissues, the tip of the ReSolve Locking Drainage Catheter is positioned within the pigtail curve. Depth markers located on the catheter shaft aid accurate placement and act as a reference for correct positioning.



Features and Benefits

Accurate placement

Depth markers aid in accurate placement and act as a reference for correct positioning.

Atraumatic

Catheter tip positioned within pigtail curve to reduce tissue irritation.

Ease of Insertion

Tapered tip and hydrophilic coating facilitate a smooth introduction to reduce patient discomfort.

Catheter Security

The unique design makes it easy to engage the locking mechanism with one hand and is also tamper resistant. The low profile of the locking mechanism and hub are designed to provide patient comfort.

Drainage Efficiency

The large drainage holes maximize drainage while maintaining catheter patency.

High Visibility

Radiopaque platinum marker band incorporated into the catheter wall for visualizing the most proximal drainage hole.

Easy Identification

Color-coded catheter hubs and packaging make it easy to determine the French size for initial placements and catheter exchanges.

Ordering Information

Catalog Number	French Size	Working Length	Marker Band	Hub Color	Drainage Holes
RLC15-6-038MB	6.5F	15 cm	Yes	Green	8
RLC15-7-038MB	7.5F	15 cm	Yes	Red	8
RLC15-8-038MB	8.5F	15 cm	Yes	Blue	8
RLC15-10-038MB	10F	15 cm	Yes	Purple	8
RLC15-12-038MB	12F	15 cm	Yes	Yellow	9
RLC15-14-038MB	14F	15 cm	Yes	Orange	9
RLC-6-038MB	6.5F	25 cm	Yes	Green	8
RLC-7-038MB	7.5F	25 cm	Yes	Red	8
RLC-8-038MB	8.5F	25 cm	Yes	Blue	8
RLC-10-038MB	10F	25 cm	Yes	Purple	8
RLC-12-038MB	12F	25 cm	Yes	Yellow	9
RLC-14-038MB	14F	25 cm	Yes	Orange	9
RLC40-6-038MB	6.5F	40 cm	Yes	Green	8
RLC40-7-038MB	7.5F	40 cm	Yes	Red	8
RLC40-8-038MB	8.5F	40 cm	Yes	Blue	8
RLC40-10-038MB	10F	40 cm	Yes	Purple	8
RLC40-12-038MB	12F	40 cm	Yes	Yellow	9
RLC40-14-038MB	14F	40 cm	Yes	Orange	9
RLC-6-038	6.5F	25 cm	-	Green	8
RLC-7-038	7.5F	25 cm	-	Red	8
RLC-8-038	8.5F	25 cm	-	Blue	8
RLC-10-038	10F	25 cm	-	Purple	8
RLC-12-038	12F	25 cm	-	Yellow	9
RLC-14-038	14F	25 cm	-	Orange	9

Packaged 5 per box

Contact Merit Medical for more information on our line of drainage products.



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
www.merit.com

Merit Medical Europe, Middle East,
and Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country
Austria 0800 295 374 | Belgium 0800 72 906 (Dutch)
0800 73 172 (French) | Denmark 80 88 00 24
France 0800 91 60 30 | Finland 0800 770 586 |
Germany 0800 182 0871 | Ireland (Republic) 1800 553 163
Italy 800 897 005 | Luxembourg 8002 25 22
Netherlands 0800 022 81 84 | Norway 800 11629
Sweden 020 792 445 | UK 0800 973 115

ID 021414 402815001/B