


# Bearing nsPVA®

Embolization Particles

## CATHETER COMPATIBILITY CHART



PART NUMBER		V100EP	V200EP	V300EP	V400EP	V600EP	V800EP	V1100EP
PARTICLE SIZE		45-150 µm	150-250 µm	250-355 µm	355-500 µm	500-710 µm	710-1000 µm	1000-1180 µm
CATHETER	INNER DIAMETER							
Angiographic Catheter	0.046" (1.17 mm)							
Angiographic Catheter	0.040" (1.02 mm)							
EmboCath® Plus Microcatheter	0.028" (0.71 mm)							
Merit Maestro® Microcatheter 2.9F	0.027" (0.69 mm)							
Merit Maestro Microcatheter 2.8F	0.024" (0.62 mm)							
SwiftNINJA® Steerable Microcatheter 2.4F	0.021" (0.54 mm)							
Merit Maestro Microcatheter 2.4F	0.020" (0.52 mm)							
Merit Pursue™ Microcatheter 2.0F	0.020" (0.50 mm)							
Merit Maestro Microcatheter 2.1F	0.018" (0.46 mm)							
Microcatheter	0.017" (0.43 mm)							
Merit Pursue Microcatheter 1.7F	0.016" (0.40 mm)							

Above catheter compatibility information based on internal testing.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

merit.com

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Singapore  
+65 67536311

South Korea  
+8210 2205 0707

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115