



## SCHEMA DI COMPATIBILITÀ DEI CATETERI

CATETERE	DIAMETRO INTERNO
Catetere angiografico	0,046" (1168 µm)
Catetere angiografico	0,040" (1016 µm)
EmboCath® Plus	0,028" (711 µm)
Merit Maestro®	0,027" (686 µm)
Merit Maestro®	0,024" (610 µm)
Merit Maestro®	0,020" (508 µm)

Dimensioni particellari compatibili



www.merit.com

Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • Paesi Bassi • Tel: +31 43 358 82 22

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Irlanda • +353 (0) 91 703 733

Numero verde per chiamate dai seguenti paesi: Austria 0800 295 374 • Belgio 0800 72 906 (olandese) 0800 73 172 (francese) • Danimarca 80 88 00 24 • Francia 0800 91 60 30

Finlandia 0800 770 586 • Germania 0800 182 0871 • Irlanda (repubblica) 1800 553 163 • Italia 800 897 005 • Lussemburgo 8002 25 22 • Paesi Bassi 0800 022 81 84 • Norvegia 800 11629

Svezia 020 792 445 • Regno Unito 0800 973 115

©2013 Merit Medical Systems, Inc. Tutti i diritti riservati.

402840004/A ID 052813