

MADE TO BREAK



FastBreak™
BREAKAWAY CONNECTOR 

FastBreak™

BREAKAWAY CONNECTOR

The FastBreak Breakaway Connector is designed to promote patient safety by having the catheter detach at the connector site rather than causing injury to the organ, soft tissues, or compromising the drainage treatment being offered.



The FastBreak Breakaway Connector is a patient safety device intended to serve as a conduit between a general drainage catheter and an external fluid reservoir. The two-piece connector has a patient side that connects to the catheter tubing and a second side that connects to a drainage supply receptacle. In use, the two sides are connected, thereby allowing fluids to flow through uninterrupted. When a force is applied to the drainage bag the connector will disconnect eliminating risk of the catheter dislodging from the patient. The connector is then able to be reconnected.

FEATURES & BENEFITS:

- Patient safety
- Prevention of bodily fluids becoming exposed or catheter being compromised
- May lessen additional services required in health care facilities
- Easy to reconnect



ORDERING INFORMATION

Part Number	Qty per Box
FBG10	10

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733