



EXCEPTIONAL INFLATION CONTROL



Q50[®] PLUS

Stent Graft Balloon Catheter

The aortic modeling balloon that offers the broadest diameter range with the lowest available profile.

Q50[®] PLUS

Stent Graft Balloon Catheter

Versatile balloon catheter ideal for abdominal and thoracic EVAR procedures. In addition to reliable AAA & TAA stent graft modeling and expansion, the Q50[®] PLUS Stent Graft Balloon Catheter is ideal for temporary occlusion of large vessels.

PRODUCT FEATURES

- Broadest balloon diameter range with the lowest available profile
 - 10mm to 50mm inflation diameters
 - 12 French introducer sheath compatible
- Strong yet compliant balloon utilizing a proprietary polyurethane technology
- Two versatile catheter lengths
 - 65cm length ideal for AAA cases
 - 100cm length ideal for TAA cases and general occlusion
- Unique tapered flexible tip facilitates device passage
- Exceptional balloon inflation and deflation rates
- Optimized catheter pushability
- Accurately positioned markerbands at balloon ends


BALLOON COMPLIANCE CHART

Recommended Inflation Volume (cc)	Balloon Diameter (mm)
3	10
6	20
16	30
32	40
44	45
60	50

ORDERING INFORMATION

Part Number	Balloon Diameter Range (mm)	Catheter Length (cm)
Q50-65P	10-50	65
Q50-100P	10-50	100

Manufactured by



Q50[®] is a registered trademark of QXMédical, LLC. QXMedical.com

Disclaimer: Products listed may not be available in all markets. As with all medical devices, please refer to the Instructions for Use. Prescription only.

Caution: USA Federal Law restricts the sale, distribution, or use of these devices to, by, or on the order of a physician.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115