



TRACKABILITY & CONFIDENCE



UAC²
UTERINE ARTERY
CATHETER


For ease and confidence
when accessing uterine arteries

NEW IMPRESS[®] UTERINE ARTERY CATHETER

The newly improved Impress Uterine Artery Catheter is designed to catheterize even the most challenging uterine vasculature.



Legato hydrophilic coating is designed to reduce friction for smoother navigation through tortuous anatomy.



Proprietary Shaft Material

for flexibility and torque without compromising steerability and kink-ability.



Catheter Hub
minimizes kinking and enhances steering.

Learn more at merit.com/UAC2.

ORDERING INFORMATION

Part Number	French Size	Length (cm)	Max Wire Size	Tip Shape	Hydrophilic Coating	Side Ports	Distal Taper	Inner Diameter	Flow Rate
59035IUACTMB-H	5	90	0.035" (0.89 mm)	UAC 2	Distal 40 cm	None	19 cm	0.046" (1.17 mm) shaft. Tip tapers to 0.035" (0.089 mm)	15 ml/sec @ 1200 PSI

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.




Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733



Extended braiding to distal portion of catheter shaft for additional support and pushability

Hydrophilic coating distal 40 cm of the catheter shaft including tip.

Tapered to 4F for improved placement when seeking vessel ostium.

New radiopaque marker band for improved visibility