




Marquis® Series Stopcocks



CONTROL • CONVENIENCE • QUALITY

Marquis® Series Stopcocks



The Marquis® Series Stopcock is another innovative accessory product from Merit. Our larger, easy-grip handle makes a Marquis Stopcock effortless to grip and turn, and our patented airless rotator prevents air bubbles.

Using the same state-of-the-art design that makes Merit's manifold the industry standard, the Marquis Series Stopcock offers unsurpassed feel and ease-of-use. Available in pressure ratings of 200, 500 and 1050 PSI with easy-to-turn handles.

Catalog Number	Body Style	Handle Color	Handle Type	Pressure Rating (psi/bar)	Luer
L3SNC	3-Way	White	Off	50 psi (3.5 bar)	Swivel Nut
L4SNC	4-Way	White	Off	50 psi (3.5 bar)	Swivel Nut
S1FLC	1-Way	White	On	200 psi (14 bar)	Fixed Male Luer
S1SLC	1-Way	White	On	200 psi (14 bar)	Slip Luer
S1SNC	1-Way	White	On	200 psi (14 bar)	Swivel Nut
S3FC	3-Way	White	Off	200 psi (14 bar)	All Female
S3FLC	3-Way	White	Off	200 psi (14 bar)	Fixed Male Luer
S3SLC	3-Way	White	Off	200 psi (14 bar)	Slip Luer
S3SNC	3-Way	White	Off	200 psi (14 bar)	Swivel Nut
S4FLC	4-Way	White	Off	200 psi (14 bar)	Fixed Male Luer
S4SLC	4-Way	White	Off	200 psi (14 bar)	Slip Luer
S4SNC	4-Way	White	Off	200 psi (14 bar)	Swivel Nut
M1SNC	1-Way	Light Blue	On	500 psi (34 bar)	Swivel Nut
M3FLC	3-Way	Light Blue	Off	500 psi (34 bar)	Fixed Male Luer
M3RLC	3-Way	Light Blue	Off	500 psi (34 bar)	Left Rotating Male Luer
M3RRC	3-Way	Light Blue	Off	500 psi (34 bar)	Right Rotating Male Luer
M3RRC-R	3-Way	Red	Off	500 psi (34 bar)	Right Rotating Male Luer
M3SLC	3-Way	Light Blue	Off	500 psi (34 bar)	Slip Luer
M3SNC	3-Way	Light Blue	Off	500 psi (34 bar)	Swivel Nut
M4FLC	4-Way	Light Blue	Off	500 psi (34 bar)	Fixed Male Luer
M4RRC	4-Way	Light Blue	Off	500 psi (34 bar)	Right Rotating Male Luer
M4SNC	4-Way	Light Blue	Off	500 psi (34 bar)	Swivel Nut
M1FLC	1-Way	Light Blue	On	500 psi (34 bar)	Fixed Male Luer
M1SLC	1-Way	Light Blue	On	500 psi (34 bar)	Slip Luer
H1FLC	1-Way	Dark Blue	On	1050 psi (72 bar)	Fixed Male Luer
H1RC	1-Way	Dark Blue	On	1050 psi (72 bar)	Right Rotating Male Luer
H1SLC	1-Way	Dark Blue	On	1050 psi (72 bar)	Slip Luer
H1SNC	1-Way	Dark Blue	On	1050 psi (72 bar)	Swivel Nut
H3FLC	3-Way	Dark Blue	Off	1050 psi (72 bar)	Fixed Male Luer
H3RLC	3-Way	Dark Blue	Off	1050 psi (72 bar)	Left Rotating Male Luer
H3RRC	3-Way	Dark Blue	Off	1050 psi (72 bar)	Right Rotating Male Luer
H3SLC	3-Way	Dark Blue	Off	1050 psi (72 bar)	Slip Luer
H3SNC	3-Way	Dark Blue	Off	1050 psi (72 bar)	Swivel Nut
H4FLC	4-Way	Dark Blue	Off	1050 psi (72 bar)	Fixed Male Luer
H4RRC	4-Way	Dark Blue	Off	1050 psi (72 bar)	Right Rotating Male Luer
H4SNC	4-Way	Dark Blue	Off	1050 psi (72 bar)	Swivel Nut
U3RRC	3-Way	Red	Off	1200 psi (83 bar)	Right Rotating Male Luer



S1SNC



M1SLC



H1FLC



U3RRC

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115