


# Mistique™

*Merit expands its line of therapeutic infusion catheters.*



**5 FRENCH MISTIQUE™ INFUSION CATHETER**

Merit provides the clinician  
cost effective infusion options.

# Mistique™

*A new 5 French infusion catheter with uniform distribution of therapeutic solutions from proximal to distal end of the infusion segment.*

## CATHETER

Manufactured to provide the proper stiffness for good trackability and ease of placement. Side holes rotate around the catheter to create a more uniform dispersion of the therapeutic agent.




### TWO TIP SIZES

The Mistique Infusion Catheter is available in two configurations that accept an 0.035" (0.89mm) wire or an 0.038" (0.97mm) wire.



### ACCESSORIES

Accessories include the Squirt™ Fluid Dispensing System and/or Pulse Infusion Kits.



### COST EFFECTIVE

The Mistique Infusion Catheter is designed to be used with either a standard guide wire or an infusion wire. Therapeutic solution is delivered uniformly throughout the entire infusion segment.

## ORDERING INFORMATION

Catalog Number	French Size	Catheter Length (cm)	Infusion Length (cm)	Recommended Guide Wire
MIC5-45-5-035	5F	45	5	0.035" (0.89mm)
MIC5-45-10-035	5F	45	10	0.035" (0.89mm)
MIC5-45-20-035	5F	45	20	0.035" (0.89mm)
MIC5-90-5-035	5F	90	5	0.035" (0.89mm)
MIC5-90-10-035	5F	90	10	0.035" (0.89mm)
MIC5-90-20-035	5F	90	20	0.035" (0.89mm)
MIC5-135-5-035	5F	135	5	0.035" (0.89mm)
MIC5-135-10-035	5F	135	10	0.035" (0.89mm)
MIC5-135-20-035	5F	135	20	0.035" (0.89mm)
MIC5-45-5-038	5F	45	5	0.038" (0.97mm)
MIC5-45-10-038	5F	45	10	0.038" (0.97mm)
MIC5-45-20-038	5F	45	20	0.038" (0.97mm)
MIC5-90-5-038	5F	90	5	0.038" (0.97mm)
MIC5-90-10-038	5F	90	10	0.038" (0.97mm)
MIC5-90-20-038	5F	90	20	0.038" (0.97mm)
MIC5-135-5-038	5F	135	5	0.038" (0.97mm)
MIC5-135-10-038	5F	135	10	0.038" (0.97mm)
MIC5-135-20-038	5F	135	20	0.038" (0.97mm)

Packaged in pouch; 5 per box.

Eto Sterilized

## OPTIONAL ACCESSORIES

K12-00462	Pulse Infusion Kit consists of in-line check relief valve, AccessPLUS™ hemostasis valve, 1ml infusion syringe, and 20ml reservoir syringe.
K12-00463	Pulse Infusion Kit consists of in-line check relief valve, 1ml infusion syringe, and 20ml reservoir syringe.
FDS-100	Squirt™ Fluid Dispensing System
MAP150	AccessPLUS™ Hemostasis Valve

**MERITMEDICAL**  
Innovations in Cardiology & Radiology

Merit Medical Systems, Inc. • 1600 West Merit Parkway, • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT  
Merit Medical International • 136, Avenue Joseph Kessel • 78960 Voisins-le-Bretonneux-France • 33 (0) 1 30 57 24 24  
Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • 353 (0) 91 703 733  
Free phone from specific country:

www.merit.com

**Belgium** 0800 72906 • **France** 0800 91 60 30 • **Germany** 0800 182 0871 • **Luxembourg** 0800 22522 • **Netherlands** 0800 022 8184 • **UK** 0800 973 115  
U.S. Patent Numbers 5,968,017 5,957,901 5,651,170 5,616,203 and 5,456,676 Other U.S. and Foreign Patents Pending 400410001/A ID 030200