



The Right Combination

*Everything necessary for the procedure is on
your shelf in an organized, ready-to-use tray*

PERICARDIOCENTESIS KIT™

BE PREPARED!

THE RIGHT COMBINATION

Everything necessary for the procedure is on your shelf in an organized, ready-to-use tray.

TIME SAVINGS

Time previously spent assembling a kit from scratch can now be focused toward the patient.

COST SAVINGS

Ordering the ready-to-use kit saves on personnel time and lowers overall component cost.

PROCEDURE SPECIFIC DESIGN

High-flow catheters are developed with pericardiocentesis procedures in mind.



8.3F Kit COMPONENTS:

- 8.3F (40 cm) pigtail or straight catheter (with marker band)
- 8.5F (22 cm) dilator
- .035" (0.89 mm) x 80 cm, double ended (J & straight) standard PTFE guide wire
- 25G x 1 1/2" (3.8 cm) anesthesia needle
- 18G x 1 1/2" (3.8 cm) medication needle
- 18G x 3.5" (9 cm) access needle
- 18G x 5.9" (15 cm) access needle
- #11 Futura® safety scalpel
- 10 mL luer lock syringe
- 20 mL luer lock syringe
- 60 mL luer lock syringe
- 1000 mL disposal bag with drainage spout
- (2) Large bore three-way stopcocks
- Double-ended alligator clip cable
- (10) 4" x 4" (10 cm x 10 cm) gauze sponges
- 60 cm x 60 cm fenestrated drape
- 46 cm x 60 cm absorbent barrier tri-layer drape
- Suture 3-0 silk with curved cutting needle
- PAL™ label set – (4) yellow saline, (4) red specimen, (4) blue lidocaine
- 2 ChloroPrep®* 3 mL (orange tint) for use in PC801 and PC802



6F Kit COMPONENTS:

- 6F (60 cm) high-flow pigtail or straight catheter
- 6F (22 cm) dilator
- .035" (0.89 mm) x 80 cm, double ended (J and straight) heavy duty PTFE guidewire
- 22G x 2 3/8" (6 cm) anesthesia needle
- 19G x 1 1/2" (3.8 cm) medication needle
- 18G x 3.5" (9 cm) access needle
- 18G x 5.9" (15 cm) access needle
- #11 scalpel
- 10 mL luer lock syringe
- 60 mL luer lock syringe
- 1400 mL disposal bag with connecting tube
- Three-way stopcock
- Double-ended alligator clip cable
- (10) 4" x 4" (10 cm x 10 cm) gauze sponges
- 40 cm x 40 cm fenestrated drape
- 46 cm x 60 cm absorbent barrier tri-layer drape

* ChloroPrep® is a registered trademark of CareFusion 2200, Inc.

Cat. No.	Description
----------	-------------

PC801	8.3F Pigtail Catheter
PC802	8.3F Straight Catheter
PC803EU	8.3F Pigtail Catheter, Europe
PC804EU	8.3F Straight Catheter, Europe
PC101	6F Pigtail Catheter
PC201	6F Straight Catheter
PC101EU	6F Pigtail Catheter, Europe
PC201EU	6F Straight Catheter, Europe

