


PERICARDIOCENTESIS TRAY



Emergency Controlled

Everything that the intervention requires is brought together in a ready-to-use tray.



PERICARDIOCENTESIS TRAY

BE PREPARED!

Merit's Pericardiocentesis trays offer:

THE RIGHT COMBINATION

Everything necessary for the procedure is on your shelf in an organized, ready-to-use tray.

TIME SAVINGS

Time previously spent assembling a kit from scratch can now be focused toward the patient.

COST SAVINGS

Ordering the ready-to-use kit saves on personnel time and lowers overall component cost.

PROCEDURE SPECIFIC DESIGN

High-flow catheters are developed with pericardiocentesis procedures in mind.

8.3F KIT COMPONENTS:

8.3F (40 cm) pigtail or straight catheter (with marker band)
8.5F (22 cm) dilator
.035" (0.89 mm) x 80 cm,
double ended (J & straight)
standard PTFE guide wire
25G x 1 1/2" (3.8 cm) anesthesia needle
18G x 1 1/2" (3.8 cm) medication needle
18G x 3.5" (9 cm) access needle
18G x 5.9" (15 cm) access needle
#11 Futura® safety scalpel
10 mL luer lock syringe
20 mL luer lock syringe
60 mL luer lock syringe
1000 mL disposal bag with drainage spout
(2) Large bore three-way stopcocks
Double-ended alligator clip cable
(10) 4" x 4" (10 cm x 10 cm) gauze sponges
60 cm x 60 cm fenestrated drape
46 cm x 66 cm absorbent barrier tri-layer drape
Suture 3-0 silk with curved cutting needle
PAL™ label set –
(4) yellow saline, (4) red specimen, (4) blue lidocaine



6F KIT COMPONENTS:

6F (60 cm) high-flow pigtail or straight catheter
6F (22 cm) dilator
.035" (0.89 mm) x 80 cm,
double ended (J and straight)
heavy duty PTFE guidewire
22G x 2 3/8" (6 cm) anesthesia needle
18G x 1 1/2" (3.8 cm) medication needle
18G x 3.5" (9 cm) access needle
18G x 5.9" (15 cm) access needle
#11 scalpel
10 mL luer lock syringe
60 mL luer lock syringe
1400 mL disposal bag with connecting tube
Three-way stopcock
Double-ended alligator clip cable
(10) 4" x 4" (10 cm x 10 cm) gauze sponges
61 cm x 61 cm fenestrated drape
46 cm x 66 cm absorbent barrier tri-layer drape



ORDERING INFORMATION

Catalogue Number	Description
PC803EU	8.3F Pigtail Catheter
PC804EU	8.3F Straight Catheter
PC101	6F Pigtail Catheter
PC201	6F Straight Catheter



Marker Band

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733