



A GLOBAL LEADER IN INFLATION

PRECISION. VERSATILITY. EASABILITY.



Inflation Devices

INFLATION DEVICES

For nearly 30 years, Merit Medical has been a market leader in inflation device technology. By working closely with clinicians, we are able to create innovative devices to meet the need for faster inflation, greater capacity, accuracy, and reliability.

basixTOUCH™

Inflation Device

FAST AND VERSATILE – THE BASIXTOUCH™

provides expanded capacity for inflation. With 35 ATM pressure capability and 30 mL syringe volume, the basixTOUCH is ideal for any interventional lab using high-pressure, large-volume balloons. With wider threads and quick-release handle, the basixTOUCH offers rapid inflation and deflation to busy labs that need top efficiency.



Catalog Number	Product Description	Units Per Box
IN8130	basixTOUCH	5
IN8152	basixTOUCH with AccessPLUS	5
IN8302	basixTOUCH with Honor	5
IN8352	basixTOUCH with Honor + Extension Tubing	5
IN8952	basixTOUCH with FLO50, insertion tool & torque device	5
IN8802	basixTOUCH with PhD	5
IN8852	basixTOUCH with PhD + Extension Tubing, insertion tool & torque device	5

Indications for Use: The basixTOUCH Inflation Device is used to inflate and deflate an angioplasty balloon or other interventional device and to measure the pressure within the balloon.


INFLATION DEVICES

Our broad product portfolio offers a wide range of features and price options—analogue and digital, low- and high-pressure capacity—for inflation, deflation, and measurement of pressures in a variety of balloons and stents.

basixTOUCH40™

Inflation Device

HIGHER PERFORMANCE WHEN YOU NEED IT MOST. The average procedure may require multiple inflations, making the time required for inflation and deflation critical. The basixTOUCH40™ delivers the speed you need and the performance you demand to meet the requirements of your busy lab.



Catalog Number	Product Description	Units Per Box
IN8140	basixTOUCH 40 Inflation Device	5

Indications for Use: The basixTOUCH40 Inflation Device is used to inflate and deflate an angioplasty balloon or other interventional device and to measure the pressure within the balloon.

INFLATION DEVICES



Analog Inflation Device

THE GOLD STANDARD INFLATION DEVICE FOR INTERVENTIONAL BALLOON INFLATION

With over 10 million devices used worldwide, busy interventional labs can rely on the basixCOMPAK™ for their procedures. The basixCOMPAK offers 30 ATM pressure and 20 mL volume and can be packaged with a wide variety of angioplasty kits and procedure packs.



Analog Readout

30 ATM

For additional kit configurations, please visit merit.com/compak.

Catalog Number	Product Description	Units Per Box
IN4130	basixCOMPAK	5
IN4330	basixCOMPAK with Honor	5
IN4352	basixCOMPAK with Honor + Extension Tubing	5
IN4530	basixCOMPAK with AccessPLUS	5
IN4802	basixCOMPAK with PhD	5
IN4852	basixCOMPAK with PhD + Extension Tubing	5
IN4952	basixCOMPAK with FLO50	5

Indications for Use: The basixCOMPAK Inflation Device is used to inflate and deflate an angioplasty balloon or other interventional device and to measure the pressure within the balloon.

INFLATION DEVICES

DIAMOND TOUCH™

Digital Inflation Device

THE LATEST INNOVATION PROVIDING HIGH PERFORMANCE, VERSATILITY, AND RELIABILITY. The DiamondTOUCH™ combines the visibility and accuracy of a digital device with the convenience of an ergonomic one-handed prep handle.



Catalog Number	Product Description	Units Per Box
IN9135	DiamondTOUCH™	5
IN9152	DiamondTOUCH™ with MAP152 - AccessPLUS™ hemostasis valve, metal wire introducer, torque device	5
IN9302	DiamondTOUCH™ with MAP302 - Honor® hemostasis valve, metal wire introducer, torque device	5
IN9352	DiamondTOUCH™ with MAP352 - Honor® hemostasis valve + bonded 20 cm tubing, metal wire introducer, torque device	5
IN9403	DiamondTOUCH™ with MAP403 - MBA™ hemostasis valve, metal wire introducer, torque device	5


Indications for Use: The DiamondTOUCH Inflation Device is used to inflate and deflate an angioplasty balloon or other interventional device and to measure the pressure within the balloon.

INFLATION DEVICES

BLUE DIAMOND™

Digital Inflation Device + Fluid Dispensing Syringe

WITH AN ACCURATE DIGITAL PRESSURE GAUGE and digital timer, you don't need to wonder how long an inflation was or what pressure was obtained. The Blue Diamond™ digital inflation device offers an integrated digital gauge and timer with the ability to record and recall prior inflation information. Ideal for ensuring accurate inflation pressure and for tracking procedure data.



Catalog Number	Product Description	Units Per Box
IN7130	Blue Diamond	5
IN7152	Blue Diamond with AccessPLUS	5
IN7352	Blue Diamond with Honor + Extension Tubing	5
IN7802	Blue Diamond with PhD	5
IN7852	Blue Diamond with PhD + Extension Tubing	5
IN7952	Blue Diamond with FLO50	5

Indications for Use: The Blue Diamond Inflation Syringe is used to inflate and deflate balloon angioplasty catheters or other interventional devices and to measure the pressure and time of inflation within the balloon during the procedure. It is also used to dispense fluids into the body and monitor the pressure of that fluid.

HEMOSTASIS VALVES AND ACCESSORIES

PATIENTS AND MEDICAL STAFF WANT TO REDUCE BLOOD LOSS AS MUCH AS POSSIBLE.

In response, Merit Medical has developed a full line of hemostasis valves to minimize blood loss in interventional medicine, while also supporting a wide array of procedures and devices. Merit Medical has been a market leader in the design and manufacture of hemostasis valves and has manufactured over 10 million valves to support the diverse needs of our users.



PhD™

Small bore 7.3F device simplifies procedures with an easy one-handed push-and-release mechanism, freeing your other hand to deliver interventional equipment.

Catalog Number	Description
MAP800	Standalone
MAP801	Metal insertion tool
MAP802	Metal insertion tool, torque device
MAP803	Plastic insertion tool
MAP804	Plastic insertion tool, Torque Device
MAP850	Bonded Extension Tubing
MAP851	Bonded Extension Tubing, Metal insertion tool, Torque Device
MAP852	Bonded Extension Tubing, Metal insertion tool, Torque Device
MAP853	Bonded Extension Tubing, Plastic insertion tool
MAP854	Bonded Extension Tubing, Plastic insertion tool, Torque Device



AccessPLUS™

Premier large bore 9F device with specially designed Tuohy valve that improves wire grip.

Catalog Number	Description
MAP150	Standalone
MAP151	Metal insertion tool
MAP152	Metal insertion tool, torque device
MAP153	Plastic insertion tool
MAP154	Plastic insertion tool, Torque Device



Honor®

Small bore 7F device with push-pull bleedback mechanism. The Honor has an innovative triple-slit bleedback seal and no Tuohy-style seal.

Catalog Number	Description
MAP300	Standalone
MAP301	Metal insertion tool
MAP302	Metal insertion tool, torque device
MAP350	Bonded Extension Tubing
MAP352	Plastic insertion tool, 20 cm bonded extension tubing, torque device
MAP353	Metal insertion tool, 50 cm bonded extension tubing, torque device
MAP354	Stopcock, Plastic insertion tool, Torque Device



MAP500
Torque Device



MAP550
Metal Insertion Tool



MAP600
Plastic Insertion Tool

