



# LIFE IMPROVED™

## Peritoneal Dialysis

MERIT IS COMMITTED  
TO DIALYSIS ACCESS

# Peritoneal Dialysis

**Flex-Neck ExxTended™**  
Upper Chest

**Flex-Neck ExxTended™**  
Upper Abdomen

**Flex-Neck® Classic**  
Lower Abdomen

**Flex-Neck ARC™**  
Lower Abdomen

The Flex-Neck® Adult catheter has an internal diameter of 3.5 mm, allowing up to 30% higher flow rates than competitors.\*

Our Flex-Neck Classic and ARC shapes are designed to accommodate various tunneling techniques and exit site options.

Infant Catheter

Pediatric Catheter

ExxTended Catheter

Flex-Neck® Straight

Flex-Neck® Classic

Flex-Neck ARC™


Catheters are available to suit multiple patient types: Obese, Adult, Adolescent, Pediatric and Infants.

\*with 2.6 mm ID

# ExxTended™ Catheter Kits


## CF-5560 AND CF-5570 COMPLETE EXXTENDED CATHETER KITS

Complete ExxTended™ Catheter Kit available in two configurations for upper abdomen and pre-sternal procedures. Unique site-to-site measuring tools designed for precise placement.



## PEDIATRIC AND INFANT CATHETERS AND ACCESSORIES

Merit Medical is committed to dialysis access and patients of all ages and sizes. We are proud to supply and support peritoneal dialysis catheters and accessories in infant, pediatric, and adolescent sizes.



\*Also available for infant CSS-242

# Laparoscopic Catheter Implantation

Merit Implantation Kits are designed for the single incision-site implantation of a PD catheter to streamline the procedure and to promote efficiency. Tunnelors and Faller Trocars aid in a tight, clean exit site; guide wires facilitate over-the-wire catheter placement. Multiple kit options are available to support your specific implantation approach and tunnelling style.

## VP-211 LAPAROSCOPIC IMPLANTATION SYSTEM

Designed to streamline deep-cuff implantation in the muscle

Cuff Implantor™ Tool

Large Dilator

Small Dilator

Luke® Guide Assembly

16 cm

Tunnelor® Tool

Unique Luke® Guide Assembly with integrated, expandable sheath for smallest opening

**Suggested for Laparoscopic and Peritoneoscopic Placement**

VP-211 Implantation System
PD Catheter Kit
CC-2300 Titanium Catheter Connector
FT-1100 Metal Faller Trocar
CS-362 Implantation Stylette (Adult)
TE-1000 (Optional) Embedding® Tool

## VP-411 LAPAROSCOPIC IMPLANTATION SYSTEM

Cuff Implantor™ Tool

Large Dilator

Luke® Guide Assembly

Plastic Faller Trocar

Longer length for patients with a higher body mass index


**Suggested for Ultrasound Guided Access and Fluoroscopic Placement**

VP-511 Implantation System
PD Catheter Kit
CC-2300 Titanium Catheter Connector
(Optional) MAK-NV™ Introducer System
(Optional) Merit Laureate® Hydrophilic Guide Wire

# Percutaneous Catheter Implantation

## VP-511 PERCUTANEOUS IMPLANTATION SYSTEM

With Antegrade Catheter Faller Trocar



**Suggested for Ultrasound Guided Access and Fluoroscopic Placement**

VP-511 Implantation System

PD Catheter Kit

CC-2300 Titanium Catheter Connector


(Optional) MAK-NV™ Introducer System

(Optional) Merit Laureate® Hydrophilic Guide Wire

## Catheter Embedding and Implantation Accessories

### TE-1000

Embedding® Tool



### FT-1100

Metal Faller Trocar



### CC-2300

Titanium Connector (Two-Piece)



### CC-1301

Plastic connector and cap



### CE-1400

Catheter extender (repair kit)



### CSS-362

Implantation Stylette





# Ordering Information

## PD Catheters\* (includes catheter, implantation stencil, connector and cap)

Product Code	Size	Number of Cuffs	O.D. (mm)	I.D. (mm)	Distance Between Cuffs (cm)	Cuff to Coil (cm)	Catheter Length (cm)
CF-5250 Classic	Adult Small	2	5.3	3.5	5	14.0	62
CF-5260 Classic	Adult Standard	2	5.3	3.5	6	15.5	62
CF-5270 Classic	Adult Large	2	5.3	3.5	7	16.5	62
CF-5450 ARC	Adult Small	2	5.3	3.5	6	14.0	62
CF-5460 ARC	Adult Standard	2	5.3	3.5	6	15.5	62
CF-5470 ARC	Adult Large	2	5.3	3.5	6	16.5	62
CF-5161 Classic	Adult Standard Straight	2	5.3	3.5	6	15.5 Distal Cuff to Catheter Tip	52
CF-5243 Classic	Adolescent	2	5.3	3.5	4.5	9.5	62
CF-5244 Classic	Adolescent	1	5.3	3.5	n/a	9.5	62
CF-5249 Classic	Adolescent Large	1	5.3	3.5	n/a	13.5	62
CF-5443 ARC	Adolescent	2	5.3	3.5	4.5	9.5	52
CF-5242 Classic	Pediatric	2	5.3	3.5	4	6.5	52
CF-5247 Classic	Pediatric	1	5.3	3.5	n/a	6.5	52
CF-5442 ARC	Pediatric	2	5.3	3.5	4.5	6.5	52
CF-4230 Classic	Infant	2	3.7	2.5	3	4.5	42
CF-4235 Classic	Infant	1	3.7	2.5	n/a	4.5	42
CF-4430 ARC	Infant	2	3.7	2.5	3	4.5	42

## ExxTended™ Coiled Catheter Kit for Upper Abdomen and Upper Chest (Pre-sternal) Catheter Placements

Product Code	Number of Cuffs	O.D. (mm)	I.D. (mm)	Distance from Distal Cuff to Coil (cm)	Distance from Proximal Cuff to End of Catheter (cm)	Overall Length Straightened (cm)
CF-5560*	(1) Upper (1) Lower	5.3	3.5	16.0 (Lower Catheter)	23.0 (Upper Catheter)	Maximum 62.0 (Lower Catheter) 62.0 (Upper Catheter). Must be trimmed to fit patient physique and exit-site location.
CF-5570*	(2) Upper (1) Lower	5.3	3.5	16.0 (Lower Catheter)	23.0 (Upper Catheter)	

## Implantation Kits

Product Code	Product Description	Quantity
VP-211 VP-210	Implantation System for Peritoneoscopic or Laparoscopic Placement of PD Catheters (Luke Guide Assembly, Small and Large Dilators, Cuff Implantor, Retrograde Tunnelor)	1 Each
VP-411	Implantation System for Laparoscopic Placement of PD Catheters (Long Luke Guide Assembly, Large Dilator, Cuff Implantor, Plastic Faller Trocar for antigrade tunneling)	1 Each
VP-511 VP-511M	Implantation System for Percutaneous Placement of PD Catheters (Luke Guide Assembly, Small and Large Dilators, Cuff Implantor, 18F Peelaway Sheath Introducer, 038" 150 cm Guide Wire, 18G Needle, Syringe, Gauze, Clip, Plastic Faller Trocar (VP-511 only), Metal Faller Trocar (VP-511M only))	1 Each

## Accessories and Embedding

Product Code	Product Description	Quantity
CC-1301*	Catheter Connector & Cap Kit (for Merit Flex-Neck Catheters ONLY)	Box of 5
CC-1360*	Catheter Connector & Cap Kit - INFANT	Box of 5
CC-2300*	Catheter Connector, Titanium, 2-piece (Flex-Neck Adult, Adolescent and Pediatric Catheters ONLY)	1 Each
CE-1400*	Repair Kit, for Flex-Neck Adult and Infant Catheters (External)	Box of 5
CS-242	Catheter Implantation Stylette for Flex-Neck Infant - 42 cm NON-STERILE ONLY	1 Each
CS-352	Catheter Implantation Stylette for Pediatric Coiled and Adult Straight Flex-Neck Catheters - 52 cm NON-STERILE ONLY	1 Each
CS-362	Catheter Implantation Stylette for Flex-Neck Adolescent & Adult Coiled Catheters - 62 cm NON-STERILE ONLY	1 Each
CSS-242	Catheter Implantation Stylette for Flex-Neck Infant - 42 cm STERILE ONLY	1 Each
CSS-352	Catheter Implantation Stylette for Pediatric Coiled and Adult Straight Flex-Neck Catheters - 52 cm STERILE ONLY	1 Each
CSS-362	Catheter Implantation Stylette for Flex-Neck Adolescent & Adult Coiled Catheters - 62 cm STERILE ONLY	1 Each
FT-1100	Metal Faller Trocar - used to create a subcutaneous antegrade tunnel	1 Each
IA-102	Air Insufflation Set - Designed for controlled infusion of filtered air into the abdomen	Box of 5
TE-1000	Embedding Tool - Used during initial catheter implantation to embed a catheter under the skin for later use	1 Each
K12T-09571	PD Percutaneous Access Pack - Designed for use with the VP-511 Percutaneous kit - sold separately.	1 Each

\* MRI Safety Information: MR Conditional—These products have been evaluated to determine MRI compatibility and have been found to be compatible under the conditions prescribed in the Instructions for Use.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle  
East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733