

TEMNO ELITE™

Biopsy System

72%
LARGER
SAMPLES



Built for Better Biopsy

The TEMNO Elite™ Biopsy System offers the reliability you expect with new features designed to provide superior soft tissue samples and total procedure control. An added patient-safety option minimizes the risk of complications, particularly during lung biopsies.² By listening to your feedback, Merit has modernized biopsy procedures with the TEMNO Elite—a device built for better biopsy.

Versatility for Your Procedures

12G, 14G	● Breast
16G	● Liver, kidney, soft tissue masses
18G	● Liver, kidney, thyroid, soft tissue masses
20G	● Lung, thyroid, soft tissue masses

Available in 6–20 cm lengths

Adjustable Throw

(10 mm or 20 mm)
Provides precise sampling for a variety of procedural situations.

Sharp Biopsy Needle

Designed to minimize tissue insertion force.

Sample Assist™ Feature

Reduces cumbersome sample retrieval with tweezers or a scalpel!



SUPERIOR SAMPLES

Total Core™ Biopsy Technology

Average 72% larger than competitor side-notch, semi-automatic devices!



PROCEDURE CONTROL

Sample Assist

"Hands-free" feature retrieves full tissue samples with a push of a button!



PATIENT SAFETY

Valved Coaxial Introducer

Intended to minimize the risk of air entering the pleural space and fluid leakage?

Visit merit.com/TemnoElite for product demonstration and ordering information.

CREATE BIOPSY SOLUTIONS



BIOPSY PROCEDURE KITS

- K12T-04321** | Biopsy Procedure Basic Kit
- K12T-04322** | Biopsy Procedure Complete Kit

REFERENCES: 1. Internal testing with 18G. Data on File. Marketing Claims for the TEMNO Elite Device (ECN147270, ECN152582). 2. Internal Testing. Data on File. Validated with data from Valved One-Step™ Centesis Catheter. See product brochure (401273001/C, ECN98478).

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
Main: +1.801.253.1600
Customer Service: +1.801.208.4300

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733