



TRANSRADIAL

PERIPHERAL

New Longer Lengths for
Transradial Peripheral Procedures

Performa[®]
Angiographic Catheter

150 cm
LENGTH AVAILABLE


Performa[®]

Angiographic Catheter

Merit Medical's popular Performa Catheters are now available in up to 150 cm lengths for performing transradial peripheral procedures.

INNER LUMEN SIZES

- 4F .042" (1.07 mm)
- 5F .046" (1.17 mm)
- 6F .054" (1.37 mm)






TRANSRADIAL

PERIPHERAL

BERENSTEIN

PIGTAIL

FEATURES & BENEFITS




Soft bumper tip provides atraumatic positioning within the vessel ostium



Up to 150 cm catheter length for easy navigation to the target vessel via the radial artery

Winged polycarbonate hub for easy manipulation



Wire braided for torque, stability, and control

Shaft material provides thermal stability at body temperature

Large inner lumens to facilitate contrast injection

Merit has solutions for every step of your transradial peripheral procedures.
Contact your Merit representative today.



PV ULT 1



PV ULT 2



PV ULT 3



PV ULT 4




PV ULT 4.5




PV BERN




PV PIGTAIL



PV MP1



PV MP2



PV RENAL



PV COBRA

ORDERING INFORMATION

Catalog Number	French Size	Inner Diameter (inches/mm)	Catheter Length (cm)	Tip Length (cm)	Number of Side Ports	Shaft/Tip Description	Tip Shape	Max Guide Wire Size
ULTIMATE 1								
PV412538U1	4	.042"/1.07 mm	125	2.0	0	Braided / BT	ULT 1	0.038"
PV415038U1	4	.042"/1.07 mm	150	2.0	0	Braided / BT	ULT 1	0.038"
PV512538U1	5	.046"/1.17 mm	125	2.0	0	Braided / BT	ULT 1	0.038"
PV512538U1SH	5	.046"/1.17 mm	125	2.0	1	Braided / BT	ULT 1	0.038"
PV513538U1	5	.046"/1.17 mm	135	2.0	0	Braided / BT	ULT 1	0.038"
PV515038U1	5	.046"/1.17 mm	150	2.0	0	Braided / BT	ULT 1	0.038"
PV612538U1	6	.054"/1.37 mm	125	2.0	0	Braided	ULT 1	0.038"
ULTIMATE 2								
PV412538U2	4	.042"/1.07 mm	125	2.0	0	Braided / BT	ULT 2	0.038"
PV512538U2	5	.046"/1.17 mm	125	2.0	0	Braided / BT	ULT 2	0.038"
PV512538U2SH	5	.046"/1.17 mm	125	2.0	1	Braided / BT	ULT 2	0.038"
PV612538U2	6	.054"/1.37 mm	125	2.0	0	Braided	ULT 2	0.038"
ULTIMATE 3								
PV412538U3	4	.042"/1.07 mm	125	2.0	0	Braided / BT	ULT 3	0.038"
PV512538U3	5	.046"/1.17 mm	125	2.0	0	Braided / BT	ULT 3	0.038"
PV512538U3SH	5	.046"/1.17 mm	125	2.0	2	Braided / BT	ULT 3	0.038"
PV612538U3	6	.054"/1.37 mm	125	2.0	0	Braided	ULT 3	0.038"
ULTIMATE 4.0								
PV412538U4	4	.042"/1.07 mm	125	2.0	0	Braided / BT	ULT 4.0	0.038"
PV512538U4	5	.046"/1.17 mm	125	2.0	0	Braided / BT	ULT 4.0	0.038"
PV512538U4SH	5	.046"/1.17 mm	125	2.0	1	Braided / BT	ULT 4.0	0.038"
PV612538U4	6	.054"/1.37 mm	125	2.0	0	Braided	ULT 4.0	0.038"
ULTIMATE 4.5								
PV412538U45	4	.042"/1.07 mm	125	2.0	0	Braided / BT	ULT 4.5	0.038"
PV512538U45	5	.046"/1.17 mm	125	2.0	0	Braided / BT	ULT 4.5	0.038"
PV512538U45SH	5	.046"/1.17 mm	125	2.0	1	Braided / BT	ULT 4.5	0.038"
PV612538U45	6	.054"/1.37 mm	125	2.0	0	Braided	ULT 4.5	0.038"
BERENSTEIN								
PV412538BERN	4	.042"/1.07 mm	125	8.0	0	Braided / BT	BERN	0.038"
PV413538BERN	4	.042"/1.07 mm	135	8.0	0	Braided / BT	BERN	0.038"
PV415038BERN	4	.042"/1.07 mm	150	8.0	0	Braided / BT	BERN	0.038"
PV512538BERN	5	.046"/1.17 mm	125	8.0	0	Braided	BERN	0.038"
PV513538BERN	5	.046"/1.17 mm	135	8.0	0	Braided	BERN	0.038"
PV515038BERN	5	.046"/1.17 mm	150	8.0	0	Braided	BERN	0.038"
PIGTAIL								
PV411038PIG	4	.042"/1.07 mm	125	6.5	5	Braided	PIGTAIL	0.038"
PV512538PIG	5	.046"/1.17 mm	125	7.5	5	Braided	PIGTAIL	0.038"
PV612538PIG	6	.054"/1.37 mm	125	8.0	5	Non-Braided	PIGTAIL	0.038"
MULTIPURPOSE 1								
PV412538MP1	4	.042"/1.07 mm	125	5.8	0	Braided/BT	MP 1	0.038"
PV512538MP1	5	.046"/1.17 mm	125	5.8	0	Braided/BT	MP 1	0.038"
MULTIPURPOSE 2								
PV412538MP2SH	4	.042"/1.07 mm	125	5.8	2	Braided/BT	MP 2	0.038"
PV512538MP2SH	5	.046"/1.17 mm	125	5.8	2	Braided/BT	MP 2	0.038"
PV612538MP2SH	6	.054"/1.37 mm	125	5.8	2	Braided	MP 2	0.038"
RENAL DOUBLE CURVE								
PV512538REN	5	.046"/1.17 mm	125	4.0	0	Braided	RENAL	0.038"
COBRA 2								
PV512538COB2	5	.046"/1.17 mm	125	9.0	0	Braided	COBRA 2	0.038"

PV = Peripheral Vasculature; BT = Bumper Tip; MP = Multi Purpose; ULT = Ultimate

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733