

# RAD BOARD® 2

RADIAL ARM BOARD

## The Next Generation of Radial



Provides built-in radiation scatter protection



THINK RADIAL.  
DISCOVER MERIT.™



# Rad Board® 2

## Radial Arm Board

### OUR RADIAL ADVANTAGE

The Rad Board 2 is the latest addition to Merit's comprehensive radial portfolio.

- Designed with a sleek, smaller profile for clinical utility
- Radiation Scatter Protection
  - 23-25% neck height<sup>1</sup>
  - 43-44% waist height<sup>1</sup>
- Easy to use with bi-plane equipment
- Reversible and reusable
- Product dimensions: 28"/71 cm x 23"/58 cm



### RELATED PRODUCTS



#### Rad Trac™

Allows easy placement of the Rad Board 2 when patient is on the table

Fits all cath lab/radiology tables

Reusable

Product dimensions:  
26.5"/67 cm x 11"/28 cm



#### Rad Board Xtra™

Allows 90° perpendicular extension of arm during access

Provides an additional 11 inches (28 cm) of working space when placed next to Rad Board 2

Reversible and reusable

Product dimensions:  
26.5"/67 cm x 11"/28 cm



#### Rad Rest®

Provides wrist and elbow support

Stabilizes and positions arm

Easily adjusts to accommodate any upper extremity length

Soft foam for increased patient comfort

# Rad Board® 2

## Radial Arm Board



*The Rad Board and accessories set-up provides convenience, comfort and maximizes working space.*



*The Rad Trac fits under the patient and the mattress allowing for the Rad Board 2 to conveniently slide under the patient with ease.*

*The Rad Board can be used with or without the RadTrac or Rad Board Xtra.*



*Easy to position under patient*

# Rad Board® 2

Radial Arm Board

## ORDERING INFORMATION

Catalog Number	Description	Order Unit
<b>RB050</b>	Rad Board 2 arm board	Sold individually, non-sterile
<b>RBX</b>	Rad Board Xtra	Sold individually, non-sterile
<b>RT913</b>	Rad Trac	Sold individually, non-sterile
<b>RR811</b>	Rad Rest	Individually packaged, 10 units per box, non-sterile

<sup>1</sup>Data on file.



THINK RADIAL. DISCOVER MERIT.™

From set up to hemostasis, Merit Medical offers a complete radial product portfolio to help you achieve the procedural **EFFICIENCY** and **SAFETY** that is essential for delivering the best possible **CARE**.



SET-UP & ACCESS



ANGIOGRAPHY



INTERVENTION



HEMOSTASIS



For more information on our radial product portfolio, scan the QR code or visit [merit.com/tr](http://merit.com/tr). For more information regarding our extensive radial education opportunities, visit [thinkradial.com](http://thinkradial.com)

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

[Merit.com](http://Merit.com)

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52  
(Italiano)  
+41 22 518 00 35

UK  
0800 973 115