


Merit PAK™

Pedal Access Kits



Your
Choice
for Limb
Salvage

Merit PAK™ Pedal Access Kits

Pedal access is a promising alternative access site for recanalization procedures in patients with Peripheral Artery Disease (PAD) and Critical Limb Ischemia (CLI)¹. Merit Medical has designed efficient kits to help physicians achieve Pedal Access: Introducing Merit PAK Pedal Access Kits.

Accurate Access

- **Echo-enhanced needles** for visibility under fluoroscopy and ultrasound
- **Clear hub needles** for visibility of flash

Mini Pedal Access Kit



Easy Insertion

- **Seamless transitions** from dilator to introducer for smooth insertions
- **Larger 4F (ID) PreludeEASE™** hydrophilic sheath for insertion of target treatment devices
- **Stiffened dilator** configuration available for tortuous anatomy

Merit PAK™ - a suite of pedal access kits including everything you need to access the arteries below the knees. The micro-needles, introducers, and dilators provide quick and easy access for limb salvage procedures.



The image displays the components of the Hydrophilic Pedal Access Kit against a white background. On the left, a dynamic splash of water is shown. In the center and right, several medical instruments are arranged: a long, thin micro-needle with a red and white handle; a shorter introducer with a red T-shaped connector and a white cap; a dilator with a red and white handle; and a clear plastic syringe with a needle. A curved, clear plastic tube is also visible, connected to the handle of the micro-needle.

Hydrophilic Pedal Access Kit

Placement & Hemostasis


- **Super elastic Nitinol core** wire provides flexibility and kink resistance²
- **Hemostasis Valve Adapter** for management of injection, aspiration, and hemostasis



Merit PAK™ Pedal Access Kits


ORDERING INFORMATION

Catalog Number	Type	Introducer OD/ID	Sheath Working Length	Needle Gauge	Needle Length	Needle Hub & Type	Guide Wire Diameter	Guide Wire Length	Hemostasis Valve
PAKM4NPD4	Mini	4/3F	10 cm	21	4 cm	Clear/Echo	0.018"	40 cm	Yes
PAKM4NPD7	Mini	4/3F	10 cm	21	7 cm	Clear/Echo	0.018"	40 cm	Yes
PAKM5NPD4	Mini	5/3F	10 cm	21	4 cm	Clear/Echo	0.018"	40 cm	Yes
PAKM5NPD7	Mini	5/3F	10 cm	21	7 cm	Clear/Echo	0.018"	40 cm	Yes
PAKS4NPD4	Stiffened	4/3F	10 cm	21	4 cm	Clear/Echo	0.018"	40 cm	Yes
PAKS4NPD7	Stiffened	4/3F	10 cm	21	7 cm	Clear/Echo	0.018"	40 cm	Yes
PAKS5NPD4	Stiffened	5/3F	10 cm	21	4 cm	Clear/Echo	0.018"	40 cm	Yes
PAKS5NPD7	Stiffened	5/3F	10 cm	21	7 cm	Clear/Echo	0.018"	40 cm	Yes
PAKPHR47NPD4	Hydrophilic	6/4F	7 cm	21	4 cm	Clear/Echo	0.018"	40 cm	Yes
PAKPHR47NPD7	Hydrophilic	6/4F	7 cm	21	7 cm	Clear/Echo	0.018"	40 cm	Yes



MINI PEDAL ACCESS KIT

- Echo-enhanced, clear hub needle
- Introducer & dilator
- Nitinol/Platinum (NPD) Guide wire
- Hemostasis Valve Adapter



HYDROPHILIC PEDAL ACCESS KIT

- Echo-enhanced, clear hub needle
- Hydrophilic sheath introducer with hemostasis valve and a side arm
- Nitinol/Platinum (NPD) Guide wire

STIFFENED PEDAL ACCESS KIT

- Echo-enhanced, clear hub needle
- Introducer & stiffened dilator
- Nitinol/Platinum (NPD) Guide wire
- Hemostasis Valve Adapter

¹ El-Sayed H.F. Retrograde Pedal/Tibial Artery Access for Treatment of Infragenicular Arterial Occlusive Disease. Methodist Debaquey Cardiovasc J. 2013 Apr-Jun; 9(2): 73–78.

² Stockel D. Nitinol – Material with Unusual Properties. Endovascular Update, Issue 1. 1998

For a complete inventory please visit our website or talk to your Merit Medical Professional.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733