



Catheter Specifications

Part Number	Inner Diameter	Outer Diameter Proximal Portion	Outer Diameter Distal Portion	Usable Length	Hydrophilic Coating Length	Max. Guide Wire	Maximum Injection Pressure	Distance From Tip to 1st Markerband	Distance Between 1st and 2nd Markerband	Distance Between 2nd and 3rd Markerband	Maximum Microsphere Size	Maximum Particle Size	Maximum Coil Size	Priming Vol.
MIV-20500	0.021" (0.54 mm)	2.9F (0.97 mm)	2.4F (0.80 mm)	125 cm	80 cm	0.018" (0.46 mm)	6,900 kPa (1,000 psi)	0.5 mm	13.5 mm	-	≤ 700 µm	≤ 710 µm	≤ 0.018" (46 mm)	0.49 mL
26SN24125ST	0.021" (0.54 mm)	2.6F (0.87 mm)	2.4F (0.80 mm)	125 cm	80 cm	0.018" (0.46 mm)	6,900 kPa (1,000 psi)	0.5 mm	13.5 mm	12.9 mm	≤ 700 µm	≤ 710 µm	≤ 0.018" (46 mm)	0.53 mL
26SN24150ST	0.021" (0.54 mm)	2.6F (0.87 mm)	2.4F (0.80 mm)	150 cm	100 cm	0.018" (0.46 mm)	6,900 kPa (1,000 psi)	0.5 mm	13.5 mm	12.9 mm	≤ 700 µm	≤ 710 µm	≤ 0.018" (46 mm)	0.59 mL

Distal Tip Cross Section



Close-up of Steerable Tip



Microsphere Compatibility

Quadrasphere®

Microspheres

Catalog Number	Dry Size (µm)	Hydrated Size (µm)	Color Code
V125QS	20-40	80-160	Gray
V225QS	30-60	120-240	Orange
V325QS	50-100	200-400	Green

Embosphere®

Microspheres

Size Range	Color Code	Code pre-filled syringe 2 mL	Microspheres Color
40-120 µm	Orange	S120GH	-
100-300 µm	Yellow	S220GH	-
300-500 µm	Blue	S420GH	-
500-700 µm	Red	S620GH	-

EmboGold®

Microspheres

40-120 µm	Orange	S120EG	red
100-300 µm	Yellow	S220EG	red
300-500 µm	Blue	S420EG	red
500-700 µm	Red	S620EG	red

Bearing nsPVA®

Embolization Particles

Vial Size Range	Color Code	Vial Product Number
45-150 µm	Yellow	V100EP
150-250 µm	Purple	V200EP
250-355 µm	Blue	V300EP
355-500 µm	Green	V400EP
500-710 µm	Orange	V600EP



Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733