

# basixTOUCH40™

Inflation Device

Up to  
**40**  
ATM



Provides versatility for  
higher pressure procedures



Offers the convenience of  
one-handed preparation

Higher performance when you need it most.

# basixTOUCH40™

## Inflation Device

### SPEED + PERFORMANCE

The average procedure may require multiple inflations, making the time required for inflation and deflation critical. The basixTOUCH40 delivers the speed you need and the performance you demand to meet the needs of your busy lab.

- **Quick-release handle** for easy release, even from high-pressure inflations, accelerating deflation time.
- **Ergonomic handle and large threads** ensure that high pressures can be reached rapidly—up to 40 atm with as few as 3 rotations of the handle.\*
- **One-handed preparation** - Preparation made easy with unique one handed air removal. Simply press handle against solid surface to prep syringe.


### VERSATILE

Today's clinicians have to be prepared to deliver a wide variety of stents and balloons. Many procedures today require high pressures. The versatility of the basixTOUCH40 can provide you the extra performance you need.

- 40 ATM pressure capacity allows inflation of high pressure interventional balloons.
- 30 mL syringe barrel accommodates additional fluid for large balloons.

### ROBUST AND RELIABLE

Designed for high and low pressure procedures, the basixTOUCH40 provides expanded capability with the security and confidence in performance that comes from using a Merit product.



**The basixTOUCH40**—Expanding your capacity with the innovation, quality, and reliability that Merit customers have enjoyed for over 25 years.

### ORDERING INFORMATION

Part Number	Description
IN8140	BasixTOUCH 40 Inflation device

*Since 1987, Merit Medical has been designing and improving inflation devices to meet the need for faster inflation, greater capacity, versatility, and reliability.*

\*Rotations vary depending on balloon volume.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.208.4300  
1.800.35.MERIT

[Merit.com](http://Merit.com)

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 225180252  
(Français)  
+41 225948000  
(Italiano)  
+41 225180035  
UK  
0800 973 115