


EVOLUTION OF A MASTERPIECE

EXPERIENCE THE GENIUS OF MERIT'S
CORONARY CONTROL SYRINGE DESIGN



CCS SYRINGES

Choose the syringe
with power and control
equal to your skills.

 **MERITMEDICAL**
Innovations in Cardiology & Radiology

CORONARY CONTROL SYRINGES

Since 1987 Merit Medical has been manufacturing and selling a variety of coronary control syringes. With years of experience, Merit has designed and created quality control syringes with features that continue to meet ever changing clinical expectations.

Clear Barrel

Provides exceptional clarity and smooth as glass feel



Solid Plunger Body

Maintains stability and durability under pressure



Rotating and Fixed Luers

Allow for flexibility and confidence in luer connection



Safety Space

Intended to minimize the potential of air bubbles introduced into the catheter



Locking Ring Option


Allows for a vacuum or negative pressure to be maintained



Smart Tip Option

Creates a unique and smooth feel during injection





- *Smaller ID barrels require less force and reduce hand fatigue during injections.*



INJECT 8™

Inject 10ⁿ™

The evolution of the coronary control syringe has now provided two revolutionary syringe designs. The Inject8 and Inject10ⁿ allow clinicians to generate more pressure with less applied force, resulting in:



- IMPROVED VISUALIZATION WITH LESS CONTRAST
- INCREASED FLOW THROUGH SMALLER FRENCH SIZED CATHETERS
- DECREASED HAND FATIGUE



Through innovative design and quality improvements, Merit Coronary Control Syringes have evolved into the tools preferred by the masters.



Ring Plunger / Ring Grip



Palm Plunger / Wing Grip




Locking Ring



Smart Tip



Palm Plunger / Ring Grip



Ring Plunger
Ring Grip
Locking Ring

Ring Plunger
Ring Grip

Palm Plunger
Ring Grip


Palm Plunger
Wing Grip

ROTATING ADAPTER

6ml Barrel	CCS460	•	•	•
Inject8™	•	CCS880	CCSB880	CCSW880
Inject10 ^{fl} ™	•	CCX010	CCXB010	•
10ml Barrel	•	CCS800	CCSB900	CCSB910
12ml Barrel	CCS100	CCS600	CCSB700	CCSB710
12ml Locking Ring	•	•	CCS200	•
20ml Barrel	•	CCS320	•	•

SMART TIP™

10ml Barrel	•	ST800	•	•
12ml Barrel	•	ST600	•	•
20ml Barrel	•	ST320	•	•



Ring Plunger
Ring Grip
Locking Ring

Ring Plunger
Ring Grip

Palm Plunger
Ring Grip

Palm Plunger
Wing Grip

FIXED MALE CONNECTOR

6ml Barrel	CCS461	•	•	•
Inject8™	•	CCS881	CCSB881	CCSW881
Inject10 ^{fl} ™	•	CCX011	CCXB011	•
10ml Barrel	•	CCS801	CCSB901	CCSB911
12ml Barrel	CCS101	CCS601	CCSB701	CCSB711
12ml Locking Ring	•	•	CCS201	•
20ml Barrel	•	CCS321	•	•

SMART TIP™

10ml Barrel	•	ST801	•	•
12ml Barrel	•	ST601	•	•
20ml Barrel	•	ST321	•	•

Packaged 25 units per box, 4 boxes per case.



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical Europe • Horsterweg 24 • 6199 AC MAASTRICHT- AIRPORT • The Netherlands • +31 - 43 - 35 88 222

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone from specific country:

Belgium 0800 72906 • France 0800 91 60 30 • Germany 0800 182 0871 • Luxembourg 0800 22522 • Netherlands 0800 022 8184 • UK 0800 973 115

U.S. Patent Numbers 5,135,488 5,449,345 5,300,02 5,425,713 5,431,62 5,458,571 5,209,732

5,118,907 5,456,676 and 5,057,078, other U.S. and Foreign Patents Pending.

401818001/A ID 102704