


# **ALIMAXX-ES™**

Fully Covered Esophageal Stent




EXPANDING THE POSSIBILITIES



**MERITMEDICAL**  
**ENDOTEK™**

# ALIMAXX-ES™

## Fully Covered Esophageal Stent




### PROXIMAL SUTURE KNOT

*Purse-string design allows repositioning and removal immediately post-placement.*



### ANTI-MIGRATION STRUTS

*Unique feature reduces stent migration.*



### FULLY COVERED

*Polyurethane cover helps decrease tissue ingrowth. Silicone lining provides a smooth inner lumen.*

### ALIMAXX-ES™ Fully Covered Esophageal Stent

Catalog Number	Stent Size Options		Delivery Device		Packaged
	Nominal Stent Diameter (mm)	Expanded Stent Length (mm)	Working Diameter (F)	Working Length (cm)	
80129-207	12	70	16	67	1 unit
80129-208	12	100	16	64	1 unit
80129-209	12	120	16	62	1 unit
80129-210	14	70	16	67	1 unit
80129-211	14	100	16	64	1 unit
80129-212	14	120	16	62	1 unit
80129-213	16	70	16	67	1 unit
80129-214	16	100	16	64	1 unit
80129-215	16	120	16	62	1 unit
80129-201	18	70	16	67	1 unit
80129-202	18	100	16	64	1 unit
80129-203	18	120	16	62	1 unit
80129-204	22	70	16	67	1 unit
80129-205	22	100	16	64	1 unit
80129-206	22	120	16	62	1 unit

**NEW SIZES**

### MAXXWIRE™ Guide Wire

Catalog Number	Guide Wire Specifications				Packaged
	O.D. (in)	Length (cm)	Shaft	Tip Style	
80129-404	0.035	180	Super Stiff	Straight	1 unit
80129-406	0.035	260	Super Stiff	Straight	1 unit

### ACCURATE PLACEMENT

- Accurate one-handed delivery over target site
- Virtually no foreshortening or elongation
- Low profile (7.4mm) pre-loaded delivery catheter
- Engineered for increased and uniform radial force



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40

Free phone for specific country:

Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71 Ireland (republic) 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 0800 6374 8633

401870001/B ID 042710