

Brighton[®]

BIPOLAR COAGULATION PROBE



EXPANDING THE POSSIBILITIES[™]



MERTMEDICA

ENDOTEK[®]

Brighton®

BIPOLAR COAGULATION PROBE



Silver Spiral Electrode

- Provides better thermal conductivity than gold*

Hydrophilic Coated Tip

- Designed to help reduce probe adherence for more effective coagulation

350 cm Length

- Ensures compatibility with most endoscopes

Dual and Single Plugs

- Fits most generators without an adapter

Purple Catheter

- Enhances endoscopic identification of the transition from tip to catheter

* Y. S. Touloukian et al. Thermal Conductivity Metallic Elements and Alloys. 1970

Catalog Number	Catheter Diameter (F)	Length (cm)	Generator Plug Type	Packaged
BP-10F-DP-350	10	350	Dual	1 unit
BP-10F-SP-350	10	350	Single	1 unit
BP-7F-DP-350	7	350	Dual	1 unit
BP-7F-SP-350	7	350	Single	1 unit

ADDITIONAL ENDOSCOPY PRODUCTS



BIG60®
Inflation Device



Elation®
3-Stage Fixed Wire
Balloon Dilator



Elation5™
5-Stage Fixed Wire
Balloon Dilator



Elation®
3-Stage Wire Guided
Balloon Dilator



TIO™
Oral Airway,
Bite Block,
O₂ Administration

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

www.merit.com/endoscopy

401985001_002 ID 041619