


# **ALIMAXX-B<sup>®</sup>**

BILIARY STENT SYSTEM



EXPANDING THE POSSIBILITIES<sup>™</sup>


 **MERRITMEDICA**  
**ENDOTEK<sup>®</sup>**

# ALIMAXX-B<sup>®</sup>

BILIARY STENT SYSTEM

## Stabilizing Sheath

Assists in accurate deployment



## Anti-Migration Struts

Unique feature designed to reduce stent migration





## Laser Cut Nitinol Construction

Virtually no foreshortening during deployment

## Anatomy Specific

High conformability aids in stent placement within tortuous anatomies and in the Hilar region





## Advanced Clinical Solutions

Optimal space between stent strut rows allows for Y-stenting while providing luminal patency


### Deployment System

- One-handed delivery mechanism
- Added precision with better control
- Low force required for deployment

# ALIMAXX-B<sup>®</sup>

BILIARY STENT SYSTEM

## FLUOROSCOPIC VIEW



This fluoroscopic marker moves proximally during deployment to indicate deployment progress

## ENDOSCOPIC VIEW



- Clear outer catheter provides increased endoscopic visibility of proximal end of stent during deployment
- Proximal end of stent is confirmed by orange endoscopic/fluoroscopic marker

Catalog Number	Stent Diameter (mm)	Stent Length (mm)	Delivery System Length (cm)	Delivery System Diameter (F)
70129-204	8	40	190	6.5
70129-205	8	60	192	6.5
70129-206	8	80	194	6.5
70129-208	10	40	190	6.5
70129-209	10	60	192	6.5
70129-210	10	80	194	6.5

## ADDITIONAL BILIARY PRODUCTS



**Monarch<sup>®</sup>**  
Digital Inflation Device



**basixCOMPAK<sup>™</sup>**  
Analog Inflation Device



**BiliQUICK<sup>™</sup>**  
Cholangiography Rapid Refill Continuous Injection Kit



**Elation<sup>®</sup>**  
3-Stage Wire Guided DASE Indicated Balloon Dilator

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733