



# EmboGold<sup>®</sup>

Microspheres

Bringing visibility to Merit Medical's  
proven microsphere technology



# EmboGold®

## Microspheres

EmboGold® Microspheres are colored to enhance visualization during handling, while still offering the controlled, targeted embolization characteristics of Merit Medical's first microsphere product, Embosphere® Microspheres. EmboGold® Microspheres are biocompatible, hydrophilic, non-resorbable, precisely calibrated microspheres, and are indicated for the embolization of hypervascularized tumors and arteriovenous malformations.

### ORDERING INFORMATION

Size Range (µm)	Color Code	1 mL* Product Number	2 mL* Product Number
40-120	Orange	--	S120EG
100-300	Yellow	S210EG	S220EG
300-500	Blue	S410EG	S420EG
500-700	Red	S610EG	S620EG
700-900	Green	S810EG	S820EG
900-1200	Purple	S1010EG	S1020EG

\*Volume of EmboGold® Microspheres supplied in sterile saline in 20 mL syringe. Sold 5 syringes per box.



### PREDICTABLE, RELIABLE DELIVERY

EmboGold® Microspheres temporarily compress for smooth microcatheter passage. Once through the microcatheter, they return to their original spherical shape and stated diameter for predictable, reliable delivery.

### NON-AGGREGATING PROPERTIES

The hydrophilic surface and spherical shape of EmboGold® Microspheres prevent aggregation in the catheter lumen and vasculature.

### DURABLE OCCLUSION

EmboGold® Microspheres are biocompatible and nonresorbable with cell adhesion properties that enable complete and durable mechanical occlusion.

### ENHANCED VISIBILITY

EmboGold® Microspheres are colored to offer better visibility during handling for enhanced safety.

### SPECIFIC SIZE RANGES

EmboGold® Microspheres are available in six specific size ranges to achieve predictable and efficient vessel occlusion.

**Before using, refer to Instructions for Use (IFU) for indications, contraindications, warnings, precautions, and directions for use.**



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.208.4300  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

merit.com

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Korea  
+8210 2205 0707

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115