



# EmboGold™

Microspheres

Bringing visibility to Merit Medical's proven microsphere technology



MERITMEDICAL  
EmboGold™  
300-500 µm  
2 mL  
2019-02-28  
STERILE  
LOT SA20EG  
X928125  
CE 2797  
LFB532007 REV 002 05/15  
4300-2002A

# EmboGold™

## Microspheres

EmboGold™ Microspheres are colored to enhance visualisation during handling, while still offering the controlled, targeted embolisation characteristics of Merit Medical's first microsphere product, Embosphere® Microspheres. EmboGold™ Microspheres are biocompatible, hydrophilic, non-resorbable, precisely calibrated microspheres, and are indicated for the embolisation of hypervascularised tumors and arteriovenous malformations.

### ORDERING INFORMATION

Size Range (µm)	Color Code	1 mL* Product Number	2 mL* Product Number
40-120	Orange	--	S120EG
100-300	Yellow	S210EG	S220EG
300-500	Blue	S410EG	S420EG
500-700	Red	S610EG	S620EG
700-900	Green	S810EG	S820EG
900-1200	Purple	S1010EG	S1020EG

\*Volume of EmboGold™ Microspheres supplied in sterile saline in 20 mL syringe. Sold 1 syringe per box.

### PREDICTABLE, RELIABLE DELIVERY

EmboGold™ Microspheres temporarily compress for smooth microcatheter passage. Once through the microcatheter, they return to their original spherical shape and stated diameter for predictable, reliable delivery.

### NON-AGGREGATING PROPERTIES

The hydrophilic surface and spherical shape of EmboGold™ Microspheres prevent aggregation in the catheter lumen and vasculature.

### DURABLE OCCLUSION

EmboGold™ Microspheres are biocompatible and nonresorbable with cell adhesion properties that enable complete and durable mechanical occlusion.

### ENHANCED VISIBILITY

EmboGold™ Microspheres are colored to offer better visibility during handling for enhanced safety.

### SPECIFIC SIZE RANGES

EmboGold™ Microspheres are available in six specific size ranges to achieve predictable and efficient vessel occlusion.



 Biosphere Medical, S.A.  
Parc des Nations - Paris Nord 2  
383 rue de la Belle Etoile  
95700 Roissy en France  
France  2797

Before using, refer to Instructions for Use (IFU) for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

[merit.com](http://merit.com)

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle  
East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733