

# Esophageal Stents

Gastroenterology

**EndoMAXX®**  
Fully Covered Esophageal Stent



**ALIMAXX-ES™**  
Fully Covered Esophageal Stent



# Esophageal Stents

## CORE ESOPHAGEAL STENT TECHNOLOGY

### Anti-Migration Struts

Designed to reduce stent migration

### Laser Cut Design

- Virtually non-foreshortening
- Improved deployment accuracy

### Advanced Patient-Focused Design

### Accurate One-Handed Deployment

- Low profile and flexible catheter
- Pre-loaded delivery system
- Pliable and tapered tip

## **EndoMAXX®**

Fully Covered Esophageal Stent

### **Unique Features**

- Migration rates as low as 11%\*
- Distal and Proximal metal sutures




## **ALIMAXX-ES™**

Fully Covered Esophageal Stent

### **Unique Features**

- Expansion force as low as 1.6 Newtons, compared to 5.4 Newtons found in woven stents\*\*
- Small diameter sizes



\* Dua et al. GIE. Efficacy and safety of a new fully covered self-expandable non-foreshortening metal esophageal stent. 2014.

\*\* Data on file



# Esophageal Stents

## ENDOMAXX® FULLY COVERED ESOPHAGEAL STENT



Catalog Number	Stent Size Options			Delivery Device		Packaged
	Stent Diameter (mm)	Stent Length (mm)	Stent Flare Diameter (mm)	Working Diameter (F)	Working Length (cm)	

### OVER-THE-WIRE

<b>MAXX-197</b>	19	70	24	24	90	1 unit
MAXX-1910	19	100	24	24	90	1 unit
MAXX-1912	19	120	24	24	90	1 unit
MAXX-1915	19	150	24	24	90	1 unit
<b>MAXX-237</b>	23	70	28	24	90	1 unit
MAXX-2310	23	100	28	24	90	1 unit
MAXX-2312	23	120	28	24	90	1 unit
MAXX-2315	23	150	28	24	90	1 unit

**UNIQUE SIZES** are noted with **RED/BOLD** part numbers

## ALIMAXX-ES™ FULLY COVERED ESOPHAGEAL STENT



Catalog Number	Stent Size Options		Delivery Device		Packaged
	Stent Diameter (mm)	Stent Length (mm)	Working Diameter (F)	Working Length (cm)	

### OVER-THE-WIRE

<b>80129-207</b>	12	70	22	67	1 unit
<b>80129-208</b>	12	100	22	64	1 unit
<b>80129-209</b>	12	120	22	62	1 unit
<b>80129-210</b>	14	70	22	67	1 unit
<b>80129-211</b>	14	100	22	64	1 unit
<b>80129-212</b>	14	120	22	62	1 unit
<b>80129-213</b>	16	70	22	67	1 unit
80129-214	16	100	22	64	1 unit
80129-215	16	120	22	62	1 unit
<b>80129-201</b>	18	70	22	67	1 unit
80129-202	18	100	22	64	1 unit
80129-203	18	120	22	62	1 unit
<b>80129-204</b>	22	70	22	67	1 unit
80129-205	22	100	22	64	1 unit
80129-206	22	120	22	62	1 unit

## ADDITIONAL STRICTURE MANAGEMENT PRODUCTS



**BIG60®**  
60 mL inflation  
device

**Elation5™**  
5-stage dilation balloon



**Elation®**  
3-stage dilation balloon



**Elation®**  
Wire guided dilation balloon



**Maxxwire®**  
Guidewire designed  
for stent placement

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.