


Getting micro guide wires into devices can be challenging.



Making Ends Meet

BOWTIE™

Guide Wire Insertion Device



GUIDE WIRE INSERTION DEVICE

The BowTie™ guide wire insertion device is an elegant solution to the problem of inserting micro guide wires into devices. Merit Medical is pleased to offer this easy, efficient product designed for “making ends meet.”



Easy—Insert the device halfway into the BowTie.



Quick—Place the BowTie and device over the proximal end of the micro guide wire.



Efficient—Remove the BowTie after “making ends meet.”

Catalog Number: BT01

Description: BowTie Guide Wire Insertion Device (maximum guide wire $\leq 0.018''/0.46\text{mm}$)

Quantity: 20 units per box (packaged individually)

Contact Merit Medical for more information about the BowTie.



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, AFRICA AND THE MIDDLE EAST (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 22

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone for specific country: Austria 0800 295 374 • Belgium 0800 72 906 (Dutch) 0800 73 172 (French) • Denmark 80 88 00 24 • France 0800 91 60 30

Finland 0800 770 586 • Germany 0800 182 0871 • Ireland (republic) 1800 553 163 • Italy 800 897 005 • Luxembourg 8002 25 22 • Netherlands 0800 022 81 84

Norway 800 11629 • Sweden 020 792 445 • UK 0800 973 115

402498001/A ID 060412