

Dilation Portfolio

Gastroenterology




Dilation Portfolio


Elation® Core Features




Rapid Deflation



Complete Deflation



Multiple Pass Durability



Additional Features

- Accurate and repeatable dilation
- Enhanced clarity
- Radiopaque markers at the shoulders of the balloon
- Catheter securement clip
- Atraumatic tip



ELATION®

3-Stage Dilation Balloon

- Built with the premium Elation Core Features
- Most aggressive price-point
- Measurement markers on the catheter
- Fixed Wire and Wireguided configurations



ELATION 5™

5-Stage Dilation Balloon

First 5-stage dilation balloon in the market

- Built with the premium Elation Core Features
- Decrease frequency of needing 2 balloons
 - reduce procedure time
 - reduce procedure cost
- Fixed Wire and Wireguided configurations



Easy To Use Inflation Devices

- Quick release trigger deflates balloon rapidly
- Consistent design and function across all 3 inflation devices
- Completely disposable
- Single piece design
- Digital gauge with timer on **Monarch®** inflation device

BIG60 ALPHA™

60mL



- Improved ergonomics
- Easier pressure relief
- Rubberized grip pad



BIG60®
60mL



basixCOMPAK™
20mL



Monarch®
20mL

Ordering Information

NEW FIXED WIRE DILATION BALLOON 5-STAGE


Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)	Packaged
EX6	5-6-7-8-9	8	2.8	180	1 unit
EX8	7-8-9-10-11	8	2.8	180	1 unit
EX10	9-10-11-12-13	8	2.8	180	1 unit
EX12	11-12-13.5-15-16	8	2.8	180	1 unit
EX15	14-15-16.5-18-19	8	2.8	180	1 unit
EX18	17-18-19-20- 21	8	2.8	180	1 unit

NEW WIREGUIDED DILATION BALLOON 5-STAGE

Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)	Packaged
EPX6	5-6-7-8-9	5.5	2.8	180	1 unit
EPX8	7-8-9-10-11	5.5	2.8	180	1 unit
EPX10	9-10-11-12-13	5.5	2.8	180	1 unit
EPX12	11-12-13.5-15-16	5.5	2.8	180	1 unit
EPX15	14-15-16.5-18-19	5.5	2.8	180	1 unit
EPX18	17-18-19-20- 21	5.5	2.8	180	1 unit
EPCX6	5-6-7-8-9	5.5	2.8	240	1 unit
EPCX8	7-8-9-10-11	5.5	2.8	240	1 unit
EPCX10	9-10-11-12-13	5.5	2.8	240	1 unit
EPCX12	11-12-13.5-15-16	5.5	2.8	240	1 unit
EPCX15	14-15-16.5-18-19	5.5	2.8	240	1 unit
EPCX18	17-18-19-20- 21	5.5	2.8	240	1 unit

UNIQUE SIZES are noted with **red/bold** part numbers

ELATION ORDERING GUIDE



Simplify inventory management with one of our 3 SKU options.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

INFLATION DEVICES

Catalog Number	Device	Inflation Device Specifications			Packaged
		Volume (mL)	Max Pressure (ATM)	Gauge	
EN-AL6012	BIG60 ALPHA	60	12	Analog	5 units
ENDO-AN6012	BIG60	60	12	Analog	5 units
ENDO-AN2030	basixCOMPAK	20	30	Analog	5 units
ENDO-DI2030	MONARCH	20	30	Digital	5 units

FIXED WIRE DILATION BALLOON 3-STAGE

Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)	Packaged
E6	6-7-8	8	2.8	180	1 unit
E8	8-9-10	8	2.8	180	1 unit
E10	10-11-12	8	2.8	180	1 unit
E12	12-13.5-15	8	2.8	180	1 unit
E15	15-16.5-18	8	2.8	180	1 unit
E18	18-19-20	8	2.8	180	1 unit

WIREGUIDED DILATION BALLOON 3-STAGE

Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)	Packaged
EPB6	6-7-8	5.5	2.8	180	1 unit
EPB8	8-9-10	5.5	2.8	180	1 unit
EPB10	10-11-12	5.5	2.8	180	1 unit
EPB12	12-13.5-15	5.5	2.8	180	1 unit
EPB15	15-16.5-18	5.5	2.8	180	1 unit
EB18	18-19-20	5.5	2.8	180	1 unit
EPCB6	6-7-8	5.5	2.8	240	1 unit
EPCB8	8-9-10	5.5	2.8	240	1 unit
EPCB10	10-11-12	5.5	2.8	240	1 unit
EPCB12	12-13.5-15	5.5	2.8	240	1 unit
EPCB15	15-16.5-18	5.5	2.8	240	1 unit
ECB18	18-19-20	5.5	2.8	240	1 unit