



PERIPHERAL INTERVENTION

POST-IMPLANTATION CARE FOR NURSES

Flex-Neck® Peritoneal Dialysis Catheters | Life Improved™

The Flex-Neck Adult catheter has an internal diameter of 3.5 mm, allowing up to 30% higher flow rates than other catheters currently on the U.S. market.


Our Flex-Neck Classic and ARC shapes are designed to accommodate various tunneling techniques and exit-site options.

Flex-Neck ExxTended™
Upper Chest

Flex-Neck ExxTended
Upper Abdomen

Flex-Neck Classic
Lower Abdomen

Flex-Neck ARC™
Lower Abdomen



Connectors & Repair Kits Reorder Chart

Catheter Manufacturer	Radiopaque Stripe Color	Catheter Size	Replacement Connector Options (Material & Part #)	Repair Kit
Merit Medical Flex-Neck Catheters Internal Diameter 3.5 (Adult Only)	Blue	Adult, Adolescent, Pediatric, Infant	Titanium Two Part CC-2300 1 each Plastic One Part CC-1301 5 per box Plastic One Part CC-1360 (Infant Only) 5 per box	CE-1400 (For all sizes of Flex-Neck Catheters) 5 per box
All Other Catheter Brands* Internal Diameter 2.6	White	All Sizes	Titanium Two Part Plastic One Part Call Catheter Manufacturer	Call Catheter Manufacturer

Peritoneal Dialysis Product Ordering Guide PD Catheters

Product Code	Size	Number of Cuffs	I.D. (mm)	Cuff to Coil (cm)	Catheter Length	Fill Volume (cc)
CF-5250 Classic	Adult Small	2	3.5	14.0	62	6.0
CF-5260 Classic	Adult Standard	2	3.5	15.5	62	6.0
CF-5270 Classic	Adult Large	2	3.5	16.5	62	6.0
CF-5450 ARC	Adult Small	2	3.5	14.0	62	6.0
CF-5460 ARC	Adult Standard	2	3.5	15.5	62	6.0
CF-5470 ARC	Adult Large	2	3.5	16.5	62	6.0
CF-5161 Classic	Adult Standard Straight	2	3.5	15.5 Distal Cuff to Catheter Tip	52	5.0
CF-5243 Classic	Adolescent	2	3.5	9.5	62	6.0
CF-5244 Classic	Adolescent	1	3.5	9.5	62	6.0
CF-5249 Classic	Adolescent Large	1	3.5	13.5	62	6.0
CF-5443 ARC	Adolescent	2	3.5	9.5	52	6.0
CF-5242 Classic	Pediatric	2	3.5	6.5	52	5.0
CF-5247 Classic	Pediatric	1	3.5	6.5	52	5.0
CF-5442 ARC	Pediatric	2	3.5	6.5	52	5.0
CF-4230 Classic	Infant	2	2.5	4.5	42	2.1
CF-4235 Classic	Infant	1	2.5	4.5	42	2.1
CF-4430 ARC	Infant	2	2.5	4.5	42	2.1


ExxTended Coiled Catheter Kit for Upper Abdomen and Upper Chest (Pre-sternal) Catheter Placements

Product Code	Number of Cuffs	O.D. (mm)	I.D. (mm)	Distance from Distal Cuff to Coil (cm)	Overall Length Straightened (cm)	Fill Volume
CF-5560	(1) Upper (1) Lower	5.3	3.5	16.0 (Lower Catheter)	Maximum 62.0 (Lower Catheter) 62.0 (Upper Catheter).	Fill volume is variable; maximum 12.0 cc. To find fill volume of trimmed, joined catheter, multiply total length of catheter (in cm) by 0.096 to get fill volume in cc.
CF-5570	(2) Upper (1) Lower	5.3	3.5	16.0 (Lower Catheter)	Must be trimmed to fit patient during procedure.	

EXTERNAL CATHETER CARE & CLEANING

All Flex-Neck Peritoneal Dialysis Catheters are made of silicone. Cleaning agents that are compatible with silicone catheters may be acceptable for use on Flex-Neck Peritoneal Dialysis Catheters.

According to the product manufacturer, electrolytically-produced sodium hypochlorite solutions (e.g., ExSept® Plus and Alcavis®50) are recommended for exit-site, connector, and catheter cleansing of all peritoneal dialysis catheters on the market.



Catheters are available to suit multiple patient types; Obese, Adult, Adolescent, Pediatric and Infants.

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.208.4300
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22
Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Austria
0800 295 374
Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)
Denmark
80 88 00 24
Finland
0800 770 586
France
0800 91 60 30

Germany
0800 182 0871
Ireland (Republic)
1800 553 163
Italy
800 897 005
Luxembourg
8002 25 22
Netherlands
0800 022 81 84

Norway
800 11629
Portugal
308 801 034
Russia
+7 495 221 89 02
Spain
+34 911238406
Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30
(Français)
+41 22 518 02 52
(Italiano)
+41 22 518 00 35
UK
0800 973 115

Merit.com



402934001/D ID 121416