

ALIMAXX-ES™

Fully Covered Esophageal Stent

The only fully covered esophageal stent available in **12, 14, 16, 18** and 22 mm diameters allowing for multitude of patient specific solutions

Laser cut design produces virtually no stent foreshortening or elongation resulting in an accurate stent placement

Proximal nylon suture for use in repositioning or removal of the stent immediate post-placement

The +5 mm proximal flare and +3 mm distal flare provide outstanding mucosal wall contact

12 mm Diameter

14 mm Diameter

16 mm Diameter

18 mm Diameter

22 mm Diameter

Optimal radial force, combined with small diameter sizes, provide an ideal solution for patients who require stent placement near the upper esophageal sphincter

Proprietary anti-migration struts face distally and are designed to minimize migration

The first FDA cleared fully covered esophageal stent

Smartzone™ technology allows for a firm mid-body with soft flared ends intended to maintain patency while reducing healthy tissue trauma

 **MERT MEDICAL**
ENDOTEK®



Merit Medical Endotek

1600 West Merit Parkway
 South Jordan, Utah 84095
 1-801-253-1600
 1-800-35-MERIT
 endotek.merit.com

**Merit Medical Europe, Middle East,
 and Africa (EMEA)**

Amerikalaan 42, 6199 AE Maastricht-Airport
 The Netherlands
 Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West
 Galway, Ireland
 +353 (0) 91 703 733

Free phone for specific country

Austria	0800 295 374	Luxembourg	8002 25 22
Belgium	0800 72 906 (Dutch)	Netherlands	0800 022 81 84
	0800 73 172 (French)	Norway	800 11629
Denmark	80 88 00 24	Sweden	020 792 445
France	0800 91 60 30	Switzerland	800 9000 5060
Finland	0800 770 586	UK	0800 973 115
Germany	0800 182 0871		
Ireland (Republic)	1800 553 163		
Italy	800 897 005		

ALIMAXX-ES™

Stent Size Options

Delivery Device

UNIQUE SIZES	Catalog Number	Mid-Body Stent Diameter (mm)	Proximal Flare Diameter (mm)	Distal Flare Diameter (mm)	Stent Length (cm)	Delivery Device Diameter (F)	Working Length
	80129-207	12	17	15	7	22	67
80129-208	12	17	15	10	22	64	
80129-209	12	17	15	12	22	62	
80129-210	14	19	17	7	22	67	
80129-211	14	19	17	10	22	64	
80129-212	14	19	17	12	22	62	
80129-213	16	21	19	7	22	67	
80129-214	16	21	19	10	22	64	
80129-215	16	21	19	12	22	62	
80129-201	18	23	21	7	22	67	
80129-202	18	23	21	10	22	64	
80129-203	18	23	21	12	22	62	
80129-204	22	27	25	7	22	67	
80129-205	22	27	25	10	22	64	
80129-206	22	27	25	12	22	62	



EXPANDING THE POSSIBILITIES™