


BAL CONVENIENCE KIT™



EXPANDING THE POSSIBILITIES™

 **MERTMEDICA**
ENDOTEK®

BAL CONVENIENCE KIT™



Merit Endotek's BAL Convenience Kit™ contains the components needed to help ensure efficiency, consistency and safety during intubated bronchoscopy requiring bronchoalveolar lavage (BAL). This all-in-one kit is ideal for busy Intensive Care Units and practices looking to streamline and standardize patient care. At Endotek we understand the value of having exactly what you need, when you need it – ensuring the best outcome for patients.



Anti-reflux tubing prevents fluid leakage when exchanging syringes and allows manual lavage return to be collected at a safe distance from the bronchoscope



TRAP CAPtain is an easy to manage bronchial wash and BAL trap, with added securement features that prevent accidental specimen loss



Labeled blue syringes and container for lidocaine use to reduce medication error.

BRC-BAL-INT-FL	QTY
Swivel Adapter	1
Bite Block	1
4x4 Gauze	10
Silicone Instrument Lube	1
1% Lidocaine, 20 mL	1
Vial Access Device	1
Luer Lock Syringe, 20 mL, Blue "Lidocaine"	1
60 mL Cup, Blue	1
Double Stick Tape	1
Slip Tip Syringe, 10 mL, Blue "Lidocaine"	2
0.9% Sodium Chloride, 250 mL	1
Luer Lock Syringe, 60 mL	2
Anti-Reflux Tubing	1
80 mL TRAP CAPtain	1
8 oz Bowl	1
80 mL Specimen Container	3
Label Set	1
Scope Reprocessing Kit	1
Endoscope Transport Drape	1

BRC-BAL-INT-SL	QTY
Swivel Adapter	1
Bite Block	1
4x4 Gauze	10
Silicone Instrument Lube	1
1% Lidocaine, 20 mL	1
Vial Access Device	1
Luer Lock Syringe, 20 mL, Blue "Lidocaine"	1
60 mL Cup, Blue	1
Double Stick Tape	1
Slip Tip Syringe, 10 mL, Blue "Lidocaine"	2
0.9% Sodium Chloride, 250 mL	1
Slip Tip Syringe, 60 mL	2
80 mL TRAP CAPtain	1
8 oz Bowl	1
80 mL Specimen Container	3
Label Set	1
Scope Reprocessing Kit	1
Endoscope Transport Drape	1

BAL Components

STERILE

Catalog Number	Description	Packaged
TC80	80 mL TRAP CAPtain	10/box
K08-03014	Anti-Reflux Tubing	20/box

Endoscope cleaning kit and transportation drape included to aid in efficient bedside pre-cleaning and safe transport of a contaminated bronchoscope.

BAL Convenience Kit™

STERILE

Catalog Number	60 mL Syringe Tip	Collection Method	80 mL Specimen Container QTY	Packaged
BRC-BAL-INT-FL	Luer Lock	Anti-Reflux Tubing	3	6/box
BRC-BAL-INT-SL	Slip Tip	80 mL Trap	3	6/box

For additional kit options, contact your local Merit Endotek sales representative to design your own custom kit.



Merit Medical Endotek

1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
merit.com/endotek

Merit Medical Europe, Middle East, and Africa (EMEA)

Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country

Austria	0800 295 374	Luxembourg	8002 25 22
Belgium	0800 72 906 (Dutch)	Netherlands	0800 022 81 84
	0800 73 172 (French)	Norway	800 11629
Denmark	80 88 00 24	Sweden	020 792 445
France	0800 91 60 30	Switzerland	800 9000 5060
Finland	0800 770 586	UK	0800 973 115
Germany	0800 182 0871		
Ireland (Republic)	1800 553 163		
Italy	800 897 005		



403147001/B 2016-03-24