

BAL CONVENIENCE KIT™



EXPANDING THE POSSIBILITIES™

BAL CONVENIENCE KIT™



Merit Endotek's BAL Convenience Kit™ contains the components needed to help ensure efficiency, consistency and safety during intubated bronchoscopy requiring bronchoalveolar lavage (BAL). This all-in-one kit is ideal for busy Intensive Care Units and practices looking to streamline and standardize patient care. At Endotek we understand the value of having exactly what you need, when you need it – ensuring the best outcome for patients.



Anti-reflux tubing prevents fluid leakage when exchanging syringes and allows manual lavage return to be collected at a safe distance from the bronchoscope



Labeled blue syringes and container for lidocaine use to reduce medication error.

Endoscope cleaning kit and transportation drape included to aid in efficient bedside pre-cleaning and safe transport of a contaminated bronchoscope.

BAL Components (Sterile)

Catalog Number	Description	Packaged
K08-03014	Anti-Reflux Tubing	20/box

BAL Convenience Kit™ (Sterile)

Catalog Number	60 mL Syringe Tip	Collection Method	80 mL Specimen Container QTY	70 mL Specimen Container QTY	Packaged
BRC-BAL-INT-FL	Luer Lock	Anti-Reflux Tubing	2	1	6/box
BRC-BAL-INT-SL	Slip Tip	70 mL Specimen Cup	2	1	6/box

BRC-BAL-INT-FL	QTY
Swivel Adapter	1
Bite Block	1
4x4 Gauze	10
Silicone Instrument Lube	1
1% Lidocaine, 20 mL	1
Vial Access Device	1
Luer Lock Syringe, 20 mL, Blue "Lidocaine"	1
60 mL Cup, Blue	1
Double Stick Tape	1
Slip Tip Syringe, 10 mL, Blue "Lidocaine"	2
0.9% Sodium Chloride, 250 mL	1
Luer Lock Syringe, 60 mL	2
Anti-Reflux Tubing	1
70mL specimen Container	1
8 oz Bowl	1
80 mL Specimen Container	2
Label Set	1
Scope Reprocessing Kit	1
Endoscope Transport Drape	1

BRC-BAL-INT-SL	QTY
Swivel Adapter	1
Bite Block	1
4x4 Gauze	10
Silicone Instrument Lube	1
1% Lidocaine, 20 mL	1
Vial Access Device	1
Luer Lock Syringe, 20 mL, Blue "Lidocaine"	1
60 mL Cup, Blue	1
Double Stick Tape	1
Slip Tip Syringe, 10 mL, Blue "Lidocaine"	2
0.9% Sodium Chloride, 250 mL	1
Slip Tip Syringe, 60 mL	2
70mL specimen Container	1
8 oz Bowl	1
80 mL Specimen Container	2
Label Set	1
Scope Reprocessing Kit	1
Endoscope Transport Drape	1

For additional kit options, contact your local Merit Endotek sales representative to design your own custom kit.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

www.merit.com/endoscopy

403147001_003 ID 110122