


SMALL DIAMETERS



BIG SOLUTIONS

ALIMAXX-ES™

Small Diameter Fully Covered Esophageal Stent

Unique clinical options intended for patient specific treatment

- Available in 12, 14 and 16mm Diameters
- Proprietary Anti-Migration Struts
- Designed for Patient Comfort

ALIMAXX-ES™

Catalog Number	Stent Size Options				Delivery Device	
	Mid-Body Stent Diameter (mm)	Proximal Flare Diameter (mm)	Distal Flare Diameter (mm)	Stent Length (cm)	Delivery Device Diameter (F)	Working Length
80129-207	12	17	15	7	22	67
80129-208	12	17	15	10	22	64
80129-209	12	17	15	12	22	62
80129-210	14	19	17	7	22	67
80129-211	14	19	17	10	22	64
80129-212	14	19	17	12	22	62
80129-213	16	21	19	7	22	67
80129-214	16	21	19	10	22	64
80129-215	16	21	19	12	22	62



EXPANDING THE POSSIBILITIES™

Merit Medical Endotek

1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
endotek.merit.com

Merit Medical Europe, Middle East, and Africa (EMEA)

Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country

Austria	0800 295 374	Luxembourg	8002 25 22
Belgium	0800 72 906 (Dutch)	Netherlands	8000 022 81 84
	0800 73 172 (French)	Norway	800 11629
Denmark	80 88 00 24	Sweden	020 792 445
France	0800 91 60 30	Switzerland	800 9000 5060
Finland	0800 770 586	UK	8000 973 115
Germany	0800 182 0871		
Ireland (Republic)	1800 553 163		
Italy	800 897 005		

