


PERFORMA®

DIAGNOSTIC CARDIOLOGY CATHETERS



Next Generation of Radial



Proven
Performance

PERFORMA® DIAGNOSTIC CARDIOLOGY CATHETERS

Introducing Merit Medical's Performa line of Diagnostic Cardiology Catheters, designed with all the torque you need, improved shaft strength for better pushability and large inner lumens for increased flow rates. The Performa line of catheters will help you reach your goal of excellent imaging in today's cardiology environment.

NYLON MATERIAL

- Kink resistant
- Optimal surface smoothness
- Thermal stability at body temperature
- Thin wall design for improved flow rate

RADIOPACITY

- Good shaft visualization under fluoroscopy
- Tip opacity ensures accurate placement and visualization

WINGED POLYCARBONATE HUB

- Enhanced handling and control

WIRE BRAID DESIGN

- Greater torque for maximum control
- Improved longitudinal stability
- High degree of maneuverability
- Provides superior pushability

TIP DESIGNS

- Excellent tip memory
- Re-forms accurately time after time

BUMPER TIP

- Added softness to the 4F Performa catheter line for patient safety

FLAT WIRE BRAID

- Increased torque
- Greater stability
- Improved kink resistance



Proven Performer | Great Radiopacity | Improved Torque

4F PERFORMA CATHETER LINE

Accepts 0.038" (0.97 mm) guide wire for physician convenience

Bumper tip softness added for safe placement of tip into ostium

Large inner diameter 0.042" (1.07 mm) provides excellent flow rates

INNER LUMEN SIZES

- 4F 0.042" (1.07 mm)
- 5F 0.046" (1.17 mm)
- 5F 0.052" (1.32 mm) high flow
- 6F 0.054" (1.37 mm)
- 6F 0.059" (1.49 mm) high flow
- Yields higher flow rates at lower pressures
- Superior guide wire trackability

STRAIN RELIEF

- Increased strength at catheter/hub interface
- Color coded for easy identification

PERFORMA®

125CM LENGTH DIAGNOSTIC CARDIOLOGY CATHETERS

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description
----------------	-------------	-------------------	--------	---------------	-------------

JUDKINS LEFT

7521-53	5	0.046"/1.17 mm	125 cm	2.5 cm	JL3.5
7521-63	5	0.046"/1.17 mm	125 cm	2.5 cm	JL4.0
7521-73	5	0.046"/1.17 mm	125 cm	2.5 cm	JL5.0
7521-83	5	0.046"/1.17 mm	125 cm	2.5 cm	JL6.0

7501-53	6	0.054"/1.37 mm	125 cm	1.9 cm	JL3.5
7501-63	6	0.054"/1.37 mm	125 cm	1.9 cm	JL4.0
7501-73	6	0.054"/1.37 mm	125 cm	1.9 cm	JL5.0
7501-83	6	0.054"/1.37 mm	125 cm	1.9 cm	JL6.0

JUDKINS RIGHT

7523-53	5	0.046"/1.17 mm	125 cm	2.5 cm	JR3.5
7523-63	5	0.046"/1.17 mm	125 cm	2.5 cm	JR4.0
7523-73	5	0.046"/1.17 mm	125 cm	2.5 cm	JR5.0
7523-83	5	0.046"/1.17 mm	125 cm	2.5 cm	JR6.0

7503-53	6	0.054"/1.37 mm	125 cm	1.9 cm	JR3.5
7503-63	6	0.054"/1.37 mm	125 cm	1.9 cm	JR4.0
7503-73	6	0.054"/1.37 mm	125 cm	1.9 cm	JR5.0
7503-83	6	0.054"/1.37 mm	125 cm	1.9 cm	JR6.0



JL 3.5



JL 4.0



JL 5.0



JL 6.0



JR 3.5



JR 4.0



JR 5.0



JR 6.0

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description
----------------	-------------	-------------------	--------	---------------	-------------

PROGRESSIVE RIGHT CORONARY (NON TORQUE)

7512-23	6	0.054"/1.37 mm	125 cm	2.5 cm	PRC
---------	---	----------------	--------	--------	-----



PRC

For additional shapes and lengths.
Please visit Merit.com.

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description	# Side Ports	Small Pigtail Curve	Angle Length
----------------	-------------	-------------------	--------	---------------	-------------	--------------	---------------------	--------------

STRAIGHT CARDIAC PIGTAIL

7529-43	5	.046"/1.17 mm	125 cm	7.5 cm	braided pig	5		
7509-43	6	.054"/1.49 mm	125 cm	7.5 cm	braided pig	5		
7569-43	6	.059"/1.49 mm	125 cm	8.0 cm	non-braided pig	5		



145 DEGREE CARDIAC PIGTAIL

7529-53	5	.046"/1.17 mm	125 cm	7.5 cm	braided pig	5	7.5 cm	
7509-53	6	.054"/1.37 mm	125 cm	7.5 cm	braided pig	5	5.5 cm	
7569-53	6	.059"/1.49 mm	125 cm	8.0 cm	non-braided pig	5	5.5 cm	



155 DEGREE CARDIAC PIGTAIL

7529-63	5	.046"/1.17 mm	125 cm	7.5 cm	braided pig	5	7.5 cm	
7509-63	6	.054"/1.37 mm	125 cm	7.5 cm	braided pig	5	5.5 cm	
7569-63	6	.059"/1.49 mm	125 cm	8.0 cm	non-braided pig	5	5.5 cm	



Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description	# Side Ports
----------------	-------------	-------------------	--------	---------------	-------------	--------------

MULTIPURPOSE



7513-C3	6	.054"/1.37 mm	125 cm	5.8 cm	MULTI A1	-
7513-53	6	.054"/1.37 mm	125 cm	5.8 cm	MULTI A2	2
7513-F3	6	.054"/1.37 mm	125 cm	5.8 cm	MULTI B1	-
7513-63	6	.054"/1.37 mm	125 cm	5.8 cm	MULTI B2	2



RADIAL APPROACH CARDIOLOGY CATHETERS

Catalog Number	French Size	Length	Side Ports	Description
MIV™ VENTRICULOGRAM CATHETER*				
4RPS110	4	110 cm	6	braided pig
4RPL110	4	110 cm	6	braided pig
4RPS110-NB	4	110 cm	6	non-braided pig
4RPL110-NB	4	110 cm	6	non-braided pig
4RPS125	4	125 cm	6	braided pig
4RPL125	4	125 cm	6	braided pig
4RPS125-NB	4	125 cm	6	non-braided pig
4RPL125-NB	4	125 cm	6	non-braided pig
5RPS110	5	110 cm	6	braided pig
5RPL110	5	110 cm	6	braided pig
5RPS110- NB	5	110 cm	6	non-braided pig
5RPL110-NB	5	110 cm	6	non-braided pig
5RPS125	5	125 cm	6	braided pig
5RPL125	5	125 cm	6	braided pig
5RPS125-NB	5	125 cm	6	non-braided pig
5RPL125-NB	5	125 cm	6	non-braided pig
6RPS110	6	110 cm	6	braided pig
6RPL110	6	110 cm	6	braided pig
6RPS110- NB	6	110 cm	6	non-braided pig
6RPL110-NB	6	110 cm	6	non-braided pig
6RPS125	6	125 cm	6	braided pig
6RPL125	6	125 cm	6	braided pig
6RPS125-NB	6	125 cm	6	non-braided pig
6RPL125-NB	6	125 cm	6	non-braided pig

*Recommended guide wire is 0.038" (0.97 mm)



Catalog Number	French Size	Internal Diameter	Length	Length of Tip	# Side Ports	Description
----------------	-------------	-------------------	--------	---------------	--------------	-------------

PERFORMA ULTIMATE™ RADIAL CATHETERS

ULTIMATE 1

410038ULT1	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100381ULT1	4	0.042" / 1.07 mm	100 cm	2.0 cm	1	bumper tip
510038ULT1	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	bumper tip
5100381ULT1	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	bumper tip
610038ULT1	6	0.054" / 1.37 mm	100 cm	2.0 cm	-	
6100381ULT1	6	0.054" / 1.37 mm	100 cm	2.0 cm	1	
5110381ULT1	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	
6110381ULT1	6	0.054" / 1.37 mm	110 cm	2.0 cm	1	
412538ULT1	4	0.042" / 1.07 mm	125 cm	2.0 cm	-	
512538ULT1	5	0.046" / 1.17 mm	125 cm	2.0 cm	-	
5125381ULT1	5	0.046" / 1.17 mm	125 cm	2.0 cm	1	



ULT1

ULTIMATE 2

410038ULT2	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100381ULT2	4	0.042" / 1.07 mm	100 cm	2.0 cm	1	bumper tip
510038ULT2	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	bumper tip
5100381ULT2	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	bumper tip
610038ULT2	6	0.054" / 1.37 mm	100 cm	2.0 cm	-	
6100381ULT2	6	0.054" / 1.37 mm	100 cm	2.0 cm	1	
5110381ULT2	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	
6110381ULT2	6	0.054" / 1.37 mm	110 cm	2.0 cm	1	
412538ULT2	4	0.042" / 1.07 mm	125 cm	2.0 cm	-	
512538ULT2	5	0.046" / 1.17 mm	125 cm	2.0 cm	-	
5125381ULT2	5	0.046" / 1.17 mm	125 cm	2.0 cm	1	



ULT2

ULTIMATE 3

410038ULT3	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100382ULT3	4	0.042" / 1.07 mm	100 cm	2.0 cm	2	bumper tip
510038ULT3	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	bumper tip
5100382ULT3	5	0.046" / 1.17 mm	100 cm	2.0 cm	2	bumper tip
610038ULT3	6	0.054" / 1.37 mm	100 cm	2.0 cm	-	
6100382ULT3	6	0.054" / 1.37 mm	100 cm	2.0 cm	2	
5110382ULT3	5	0.046" / 1.17 mm	110 cm	2.0 cm	2	
6110382ULT3	6	0.054" / 1.37 mm	110 cm	2.0 cm	2	
512538ULT3	5	0.046" / 1.17 mm	125 cm	2.0 cm	-	



ULT3

ULTIMATE 4

(4 cm between primary & secondary curve)

410038ULT4-T40	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100381ULT4-T40	4	0.042" / 1.07 mm	100 cm	2.0 cm	1	bumper tip
510038ULT4-T40	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	bumper tip
5100381ULT4-T40	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	bumper tip
610038ULT4-T40	6	0.054" / 1.37 mm	100 cm	2.0 cm	-	
6100381ULT4-T40	6	0.054" / 1.37 mm	100 cm	2.0 cm	1	
5110381ULT4-T40	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	
6110381ULT4-T40	6	0.054" / 1.37 mm	110 cm	2.0 cm	1	

(4.5 cm between primary & secondary curve)

410038ULT4-T45	4	0.042" / 1.07 mm	100 cm	2.0 cm	-	bumper tip
4100381ULT4-T45	4	0.042" / 1.07 mm	100 cm	2.0 cm	1	bumper tip
510038ULT4-T45	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	bumper tip
5100381ULT4-T45	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	bumper tip
610038ULT4-T45	6	0.054" / 1.37 mm	100 cm	2.0 cm	-	
6100381ULT4-T45	6	0.054" / 1.37 mm	100 cm	2.0 cm	1	
5110381ULT4-T45	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	
6110381ULT4-T45	6	0.054" / 1.37 mm	110 cm	2.0 cm	1	
512538ULT4-T45	5	0.046" / 1.17 mm	125 cm	2.0 cm	-	
5125381ULT4-T45	5	0.046" / 1.17 mm	125 cm	2.0 cm	1	



ULT4

IMPRESS® MATERIAL ULTIMATE

5IMP10038ULT1	5	0.046" / 1.17 mm	100 cm	2.0 cm	-	Recommended guide wire is 0.038" (0.97 mm)
5IMP100381ULT1	5	0.046" / 1.17 mm	100 cm	2.0 cm	1	
5IMP110381ULT2	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	
5IMP110381ULT4-T40	5	0.046" / 1.17 mm	110 cm	2.0 cm	1	

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description	
PERFORMA INTERNAL MAMMARY LEFT - BRACHIAL / RADIAL CATHETERS*						
7242-10	4	0.042"/1.07 mm	100 cm	4.2 cm	ImLBR1	bumper tip
7242-20	4	0.042"/1.07 mm	100 cm	4.2 cm	ImLBR2	bumper tip
PERFORMA LEFT CORONARY - BRACHIAL / RADIAL CATHETERS*						
7240-10	4	0.042"/1.07 mm	100 cm	3.4 cm	LBR1	bumper tip
7240-20	4	0.042"/1.07 mm	100 cm	3.4 cm	LBR2	bumper tip
7240-30	4	0.042"/1.07 mm	100 cm	6.0 cm	LBR3	bumper tip
PERFORMA RIGHT CORONARY - BRACHIAL / RADIAL CATHETER*						
7241-10	4	0.042"/1.07 mm	100 cm	3.4 cm	RBR	bumper tip
PERFORMA AMPLATZ LEFT CORONARY - BRACHIAL / RADIAL CATHETER*						
7703-B0	4	0.042"/1.07 mm	100 cm	1.5 cm	AL2	bumper tip



ImLBR1



ImLBR2



LBR1



LBR2



LBR3




RBR



AL2

*Recommended guide wire is 0.038" (0.97 mm)

*Recommended guide wire is 0.038" (0.97 mm)

Catalog Number	French Size	Catheters 100 Cm	Side Port	Sheath Length	Sheath Guide Wire Type	
RADIAL MULTIPACK PLUS*						
RADMPP-1	5	ULT 1		11 cm	SS/SS	
RADMPP-2	5	ULT 1		11 cm	SS/SS	
RADMPP-3	5	ULT 1	1	7 cm	SS/SS	
RADMPP-4	5	ULT 1		7 cm	NT/PT	
RADMPP-5	5	ULT 1		11 cm	NT/PT	
RADMPP-6	5	ULT 1	1	7 cm	NT/PT	
RADMPP-11	5	ULT 4T-40		7 cm	NT/PT	
RADMPP-12	5	ULT 4T-40		11 cm	NT/PT	
RADMPP-7	5	JL 3.5, JR 4.0		7 cm	SS/SS	
RADMPP-8	5	JL 3.5, JR 4.0		11 cm	SS/SS	
RADMPP-9	5	JL 3.5, JR 4.0		7 cm	NT/PT	
RADMPP-10	5	JL 3.5, JR 4.0		11 cm	NT/PT	

Contains 1) InQwire® guide wire 0.035" (0.89 mm) X 180 cm/Fixed Core 3 mm J

SS = Stainless Steel NT = Nitinol PT = Platinum Tip



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Merit.com

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg 8002
25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Spain
+34 91 1238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52
(Italiano)
+41 22 518 00 35

UK
0800 973 115