



# DiamondTOUCH™

## Digital Inflation Device

Up to  
**35**  
ATM



Offers the convenience of  
one-handed preparation



High resolution digital  
gauge provides exceptional  
visibility and accuracy

The latest innovation providing high performance, versatility, and reliability.

# DiamondTOUCH™

## Digital Inflation Device

### 35 ATM, 30 ML DIGITAL INFLATION DEVICE

The DiamondTOUCH combines the visibility and accuracy of a digital device with the convenience of an ergonomic one-handed prep handle.

### DIAMONDTOUCH FEATURES

- One button functionality for simplicity and ease of use
- Large dial face with oversized time and pressure digits to enhance visibility
- Last inflation mark recalls the highest pressure of the previous inflation
- High resolution from multiple angles
- Ergonomic one-handed prep handle

### SPEED AND PERFORMANCE

The average procedure may require multiple inflations, making the time required for inflation and deflation critical. The DiamondTOUCH delivers the speed and performance needed to meet the needs of a busy lab.

- **Rapid Inflation** - Pressures can be reached rapidly - up to 35 atm with as few 3 rotations of the handle. \*
- **Quick release ergonomic handle** - Releases easily, even from high-pressure inflations.
- **One-handed preparation** - Preparation made easy with unique one-handed air removal. Simply press the handle against a solid surface.

### VERSATILE

The DiamondTOUCH can provide the extra performance needed:

- 30 mL syringe barrel accommodates additional fluid for larger balloons
- 35 atm pressure capacity allows inflation of high pressure interventional balloons

### ROBUST AND RELIABLE

Strong components provide expanded capability with the security and confidence in performance that comes from using Merit products.

### ORDERING INFORMATION

Part Number	Description
IN9135	DiamondTOUCH
IN9152	DiamondTOUCH with MAP152 – Access Plus hemostasis valve with metal wire introducer and torque device
IN9302	DiamondTOUCH with MAP302 – Honor hemostasis valve with metal wire introducer and torque device
IN9352	DiamondTOUCH with MAP352 – Honor hemostasis valve with bonded 20cm tubing, metal wire introducer and torque device
IN9403	DiamondTOUCH with MAP403 – MBA hemostasis valve with metal wire introducer and torque device

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.208.4300  
1.800.35.MERIT

merit.com

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Korea  
+8210 2205 0707

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115